



Lasso™

Rotatable Polypectomy Snare (COLD)

Instructions for use.
Read carefully prior to use.

Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Diversatek Healthcare

Corporate Headquarters

102 E. Keefe Ave.

Milwaukee, WI 53212 USA

P 800.558.6408 | 414.265.7620

sales@diversatekhc.com

DiversatekHealthcare.com

Diversatek Healthcare Lasso™ Rotatable Polypectomy Snare (COLD)*Instructions for Use.**Read carefully prior to use.***INTENDED USE**

The Lasso™ Rotatable Polypectomy Snare (COLD) is used for endoscopic removal of foreign bodies, especially for small polyps. It is used without diathermic energy. The device is supplied sterile and intended for single use only.

| <u>Part Number</u> | <u>Description</u> |
|--------------------|---|
| 1180-12 | Oval, 10 mm, Rotatable, 230 cm, Cold Only |
| 1180-13 | Oval, 15 mm, Rotatable, 230 cm, Cold Only |
| 1180-14 | Hex, 10 mm, Rotatable, 230 cm, Cold Only |
| 1180-15 | Hex, 15 mm, Rotatable, 230 cm, Cold Only |

CONTRAINDICATIONS

Contraindications include, but are not limited to, those specific to endoscopic polypectomy.

WARNINGS AND PRECAUTIONS

1. If the product package is open or damaged when received, do not use this device.
2. The snares require an endoscope with a minimum working channel of 2.8mm.
3. Do not use this device for any purpose other than the stated intended use.
4. Verify the expiration date on the package label prior to using the product. If the expiration date is lapsed, do not use.
5. The device should only be used by a trained medical professional.
6. A full understanding of the technical principles, clinical applications, risks associated with cold snaring (non-electrical) polypectomy and tissue resection is necessary before using this product.
7. Disengage the snare from the polyps if there is risk of complication.
8. Be careful when grasping the tissue to avoid grasping tissue or organs that are not suitable for retrieval.

INSTRUCTIONS FOR USE

1. Upon removing the device from the package, uncoil the device and visually inspect for kinks, bends, breaks, fraying, or other damage. If an abnormality is detected that would prohibit proper working condition, do not use.
2. Fully retract and extend the snare to confirm smooth operation of the device.
3. When the polyp is in endoscopic view, introduce the sheath and retracted snare into the endoscope accessory channel.
4. Advance the device in small increments until it is endoscopically viewed exiting the endoscope.
5. Extend the snare loop into position around the target polyp and pull the handle to resect.
6. Retrieve the polyp and prepare the specimen per institutional guidelines.
7. Upon completion of the polypectomy, retract the snare into the sheath and remove the device from the endoscope.
8. Upon completion of the procedure, dispose of this device per institutional guidelines for biohazardous medical waste.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: Perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia, or arrest. Main symptoms were abdominal pain, fever, and short intestinal cramps. Direct observation is required. Improper snare orientation or position can cause patient injury.

Lasso™ is a trademark of Diversatek Healthcare.

Diversatek™ is a trademark of Diversatek, Inc.

DESCRIPTION OF SYMBOLS USED ON LABELS



Manufactured for



Use-by date
(YYYY-MM-DD)



Lot number



Part number



Sterilized using ethylene oxide



Temperature limit



Humidity limitation



Atmospheric pressure limitation



Keep dry



Do not use if package is damaged



Keep away from sunlight



Do not re-use



Consult instructions for use



Not made with natural rubber latex