# *innerVision*<sup>®</sup>

*transillumination system*

Diversatek Healthcare InnerVision®Transillumination System	Page 2
Diversatek Healthcare InnerVision®Sistema de Transiluminación	Page 7
Système de transillumination Diversatek Healthcare InnerVision®	Page 12
Sistema di transilluminazione InnerVision® Diversatek Healthcare	Page 17
InnerVision®Transilluminationssystem von Diversatek Healthcare	Page 21
Diversatek Healthcare InnerVision®doorlichtsysteem	Page 26

**Diversatek**<sup>™</sup>  
Healthcare

Corporate Headquarters  
102 E. Keefe Ave.  
Milwaukee, WI 53212 USA  
P 800.558.6408 | 414.265.7620

[sales@diversatekhc.com](mailto:sales@diversatekhc.com)

[DiversatekHealthcare.com](http://DiversatekHealthcare.com)

**Diversatek Healthcare Innervision® Transillumination System**

Instructions for Use and Cleaning.  
Read carefully prior to use.

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.  
U.S. Patents #5,624,432; 5,718,666;  
International Patents #0 568 656; 69225583.4; 659,870;  
Additional international patents pending.

FIGURE 1. INNERVISION TRANSILLUMINATION SYSTEM

**DESCRIPTION**

The Diversatek Healthcare InnerVision Transillumination System is designed to aid in the identification of the esophagus, rectum and other structures by transillumination during laparoscopy, thoracoscopy, or open procedures. When used in the esophagus, the InnerVision can also be described as a lighted bougie.

The system consists of two (2) components:

1. InnerVision Fiber Optic Cable: Consists of a silicone elastomer fiber optic cable with a built-in extension cable, having a stainless steel tip connector at one end and a light source connector at the other end. Provided non-sterile. Reusable. Packaged one (1) cable per box.
  - Available in 40 French, 50 French, and 56 French sizes
  - Available with a Circon-ACMI/Olympus light source connector, or
  - Available with a Storz light source connector
  - Wolf/Dyonics light source adapter available for purchase. Used with a Circon-ACMI/Olympus cable.
2. InnerVision Detachable Tip: Clear, flexible silicone tip that fits on the tip of the Fiber Optic Cable. Sterile, single use only. Packaged five (5) tips per box.
  - Available in 40 French, 50 French, and 56 French sizes

InnerVision Transillumination System can be used with a Cogent light source, but an adapter must be ordered from Cogent Light.

All contents of the InnerVision Transillumination System are not made with natural rubber latex.

**DESCRIPTION OF SYMBOLS USED ON LABELS**

	Non-Sterile Device
	Caution, Consult Accompanying Documents
	Fragile, Handle with Care
	Consult Operating Instructions
	Do Not Re-Use
	Do Not Use if Package is Damaged
	Do Not Re-Sterilize
	Sterilized Using Ethylene Oxide

**INDICATIONS FOR USE**

The InnerVision Transillumination System is intended for use:

1. To transilluminate the esophagus and assist in the identification and location of the esophagus and its surrounding tissues during laparoscopic, thoracoscopic and open procedures.
2. To transilluminate the rectum and rectosigmoid and assist in the identification and location of these and the surrounding tissues during gynecological and colorectal laparoscopic and open procedures.

**CONTRAINDICATIONS**

Use of the InnerVision System is contraindicated in situations where the size or presence of the device could cause injury or contamination. The device is contraindicated if the patient has a severe stricture or obstruction which would not allow any size to be inserted. The device is contraindicated if the

patient has known or suspected intubation complications and/or difficulties. Use of the InnerVision System is contraindicated if resistance is felt during insertion of this device.

### **COMPLICATIONS**

- Inability to place the InnerVision System may occur. It is possible that several strictures and obstructions may prevent insertion of the device. It is possible that a smaller InnerVision System may be inserted more easily. This device should not be used if resistance is felt during insertion of the device.
- It has been reported that non-illuminated bougies have perforated the walls of the esophagus. It may be possible to perforate tissues with the InnerVision lighted bougie.
- Active severe inflammation may interfere with the passage of the larger devices and may require the use of device with a smaller diameter.
- Active bleeding may interfere with the transillumination of the esophagus.

### **WARNINGS**

1. The Fiber Optic Cable CANNOT be used on a patient without the correctly sized Detachable Tip properly connected, as thermal injury will occur.
2. Do NOT lubricate the connection between the Fiber Optic Cable & Detachable Tip.
3. This device is NOT RECOMMENDED for use on patients with known or suspected intubation complications and/or difficulties. Discontinue use if resistance is felt during insertion of this device. Resistance may cause the Detachable Tip to become separated from the Fiber Optic Cable.
4. When inserting and removing the InnerVision system from the patient, DO NOT ALLOW SHARP BENDS. The Detachable Tips have become disconnected while in the body due to sharp bending.
5. As the InnerVision is removed from the esophagus, the Detachable Tip must NOT be bent past 45 degrees at the connector; thermal injury would then be possible through the sidewall of the tip.

### **PRECAUTION**

The Detachable Tips are sterile and single use only. Dispose the tip after use. DO NOT REUSE. Detachable Tips can not be reused due to both the silicone elastomer absorbing the disinfectant solutions and the small openings in the tip retaining some of the disinfectant solution, which subsequently would allow the chemicals to leach out and cause tissue reaction in patients.

A damaged Detachable Tip must not be used. Using a damaged Detachable tip could result in unintentional detachment of the tip or exposure of the patient to the fiber optic bundle with the possible result of injury to the patient.

The Fiber Optic Cable is a delicate device. Meticulous care must be taken to prevent damage from occurring. The cable must be handled and cleaned carefully. Improper maintenance will not only affect the function of the Diversatek Healthcare InnerVision Transillumination System, but can also cause harm to the patient. Any damage will render the cable unusable. A damaged Fiber Optic Cable must not be used and must be replaced.

1. Avoid contact with any sharp edges or pointed objects. Any cut or puncture can expose the optical fibers.
2. Do not use any solvent except isopropyl alcohol.
3. Do not soak in any solution including disinfectants such as glutaraldehyde.
4. Do not use germicides such as bichloride of mercury or other corrosive solutions.
5. Do not clean in an ultrasonic cleaning device.
6. Do not use clamps or forceps on the Fiber Optic Cable or Detachable Tips.
7. Do not twist, stretch, kink, squeeze, scratch or shake the Fiber Optic Cable. Do not pull, tightly coil or severely bend the cable as irreversibly damage to the glass fibers will occur.
8. Take care in preventing any scratching from occurring on the face of glass fibers at the connector or adapter ends of the Fiber Optic Cable.
9. Cracking, discoloration, and broken fibers are the most obvious indications of aging and will affect the safety and effectiveness of the cable.
10. The Fiber Optic Cable cannot be used on a patient without the Detachable Tip connected, as thermal injury will occur.
11. When using the Circon-ACMI light source, half of the stainless steel connector on the Circon-ACMI Fiber Optic Cable is still visible (1 inch section). Do not force the cable into the Circon-ACMI/Olympus light source any further than 1 inch. Do not touch the visible part of the stainless steel connector on the cable, as it will cause thermal injury.
12. When removing the Fiber Optic Cable from the light source, hold the black grip to pull the cable. Pulling on any other part of the cable will result in damage to the fiber optics in the cable.
13. When inserting and removing the InnerVision System from the patient, DO NOT ALLOW SHARP BENDS. The Detachable Tips have become disconnected while in the body due to sharp bending.
14. This device is NOT RECOMMENDED for use on patients with known or suspected intubation complications and/or difficulties. Discontinue use if resistance is felt during insertion of this device. Resistance may cause the Detachable Tip to become separated from the Fiber Optic Cable.

Fiber optic breakage can result in light and function loss. Hold the tip connector of the Fiber Optic Cable towards a regular light and examine the light source connector end of the cable for broken fibers. Small black dots on the surface of the glass fibers indicate broken fibers (see Figure 2). Do not use the fiber optic light source for this test as the intense light can cause eye damage. A damaged Fiber Optic Cable must not be used and must be replaced.


FIGURE 2. FIBER OPTIC CABLE REPLACEMENT CHART

### **CLEANING & STERILIZING THE FIBER OPTIC CABLE**

1. Cleaning is recommended for new and used InnerVision Fiber Optic Cables prior to use. Clean the cable thoroughly, immediately after use, with a soft-bristle brush in lukewarm soapy water to remove contaminants. Use a non-oily cleaner or mild non-abrasive soap. Do not use synthetic detergents or oil-based soaps as these may be absorbed by the cable surface and may subsequently leach out and cause tissue reaction.
2. Rinse thoroughly in lukewarm water. Follow with a thorough rinse in distilled water.

3. Sterilize the Fiber Optic Cable using a steam autoclave or STERRAD® Sterilization System. Use the “STERRAD® 100 System Operator’s Manual” when sterilizing the Fiber Optic Cable with a STERRAD® System. Whichever sterilization method is used, the institutional protocols and instructions in the appropriate owner’s manual should be followed. Diversatek Healthcare’ validated gravity-type steam sterilization settings are 132 degrees Celsius for a cycle time of 15 minutes. Prevacuum steam sterilization settings are 132 degrees Celsius for a cycle time of 4 minutes. Cables sterilized by gravity-type or prevacuum steam sterilization machines should be individually wrapped. Diversatek Healthcare relies on the user to properly validate the sterilization method selected and equipment used at the user facility.
4. If sterilization is being done by steam autoclave, allow the Fiber Optic Cable to completely cool to room temperature (25°C). Our testing indicated that a cooling time of 20 minutes at room temperature is required. Do not immerse in cold liquids. Avoid contact with the metal tips. Do not attempt to attach a Detachable Tip until it is cooled to avoid burns to user and patient.
5. Do not use any solvent except isopropyl alcohol.
6. Do not soak in any solution including disinfectants such as glutaraldehyde.
7. Do not use germicides such as bichloride of mercury or other corrosive solutions.
8. Do not clean in an ultrasonic cleaning device.
9. Do not clean or sterilize the Detachable Tips. Detachable Tips can not reused due to both the silicone elastomer absorbing the disinfectant solutions and the small openings in the tip retaining some of the disinfectant solution, which subsequently would allow the chemicals to leach out and cause tissue reaction in patients.

#### **PRIOR TO USE**

1. Ensure that a functional 250-300 watt fiber optic light source is available that will accommodate the connector on the Fiber Optic Cable. An additional adapter will need to be used with a Wolf light source or a Cogent light source. A lower-power light source or an incorrect adapter will result in a lower level of illumination.
2. Inspect the Fiber Optic Cable for defects by the following steps. Damaged cables must be discarded. A backup cable should be available at time of use.
  - a. Examine the cable for cracks, splits, gouges, and tears.
  - b. Examine the junction between the connector tips and the cable for cracks, peeling, gaps or general deterioration and loosening.
  - c. Hold the tip connector towards a regular light and examine the light source adapter end of the Fiber Optic Cable for broken fibers. Replace the cable if too many fibers are broken (see Figure 2). Damage to the fibers will reduce the amount of light transmitted so that the bougie may not provide adequate illumination. If light transmission is decreased, replace the Fiber Optic Cable.
3. Attach the Detachable Tip to the Fiber Optic Cable by grasping the tip at the socket end and the cable at the connector tip and pushing the tip socket onto the mating connector tip of the cable (see Figure 3). DO NOT lubricate the connection between the Fiber Optic Cable & Detachable Tip. The connector on the 40 Fr. tip is different from the connector on the 50 Fr. and the 56 Fr.; however, the tips are attached in exactly the same way. The Fiber Optic Cable cannot be used on a patient without the Detachable Tip connected, as thermal injury will occur.

FIGURE 3. ATTACHMENT OF THE DETACHABLE TIP TO THE FIBER OPTIC CABLE.

4. Attach the Fiber Optic Cable to the light source by plugging the connector into the light source in the same manner as a standard light source extension cable. An adapter needs to be attached to the Circon-ACMI/Olympus cable before plugging into a Wolf/Dyonics light source or Cogent light source. When using the Circon-ACMI light source, half of the stainless steel connector on the Circon-ACMI Fiber Optic Cable is still visible (1 inch section). Do not force the cable into the Circon-ACMI/Olympus light source any further than 1 inch. Do not touch the visible part of the stainless steel connector on the cable, as it will cause thermal injury.

#### **ESOPHAGEAL TRANSILLUMINATION**

1. The patient is prepared for surgery and the surgical procedure initiated. The tipped end of the InnerVision lighted bougie is passed into the esophagus typically by an anesthesiologist. The InnerVision lighted bougie must be lubricated for passage down the esophagus. Water, sterile saline or a clear medical lubricating gel can be used to lubricate the lighted bougie.
2. The light source is turned on when appropriate, but always before passage of the bougie into the stomach.
3. When inserting and removing the InnerVision System from the patient, DO NOT ALLOW SHARP BENDS. The Detachable Tips have become disconnected while in the body due to sharp bending.
4. Discontinue use if resistance is felt during insertion of this device. Resistance may cause the Detachable Tip to become separated from the Fiber Optic Cable.
5. When the surgical procedure or use of the bougie is completed, turn the light source off and disconnect the InnerVision System from the light source. Hold the black finger grip and pull the connector away from the light source. Pulling on any other part of the cable will result in damage to the fiber optics in the cable.
6. As the InnerVision is removed from the esophagus, the Detachable Tip must not be bent past 45 degrees at the connector; thermal injury would then be possible through the sidewall of the tip.
7. To remove the Detachable Tip, grasp the tip near the proximal end of the Fiber Optic Cable. With hands close together, rotate the hand holding the tip in a peeling motion (see Figure 4). Do not pull the tip horizontally as damage to the connector on the cable and damage to the fiber optics will occur (see Figure 5). The Detachable Tip is discarded and the Fiber Optic Cable is cleaned and sterilized for future use.

FIGURE 4. REMOVE THE DETACHABLE TIP WITH A PEELING MOTION

FIGURE 5. INCORRECT REMOVAL OF DETACHABLE TIP

**OTHER TRANSILLUMINATION**

The same principles apply for the rectosigmoid or other uses. It is customary to have separate InnerVision Systems for esophageal and rectal uses.

**STORAGE**

The InnerVision Fiber Optic Cable should be stored in a cool, dry place in temperatures not exceeding 25°C. The original carton lined with foam provides an adequate storage environment for the Fiber Optic Cable. Do not tightly coil or kink the cable while storing the InnerVision System.

**RETURNED GOODS POLICY**

Authorization for the return of merchandise must be obtained from Diversatek Healthcare, Inc. Merchandise returned must have all manufacturer’s seals intact and be received within 30 days of invoice to be eligible for credit or replacement. Returned products may be subject to restocking charges.

**DIVERSATEK HEALTHCARE WARRANTY**

Limited Warranty, Limitation of Liability and Disclaimer. Diversatek Healthcare warrants to the purchaser that: (a) the InnerVision Transillumination Fiber Optic Cable will be free from defects in materials and workmanship for a period of twelve (12) months from the time of shipment. Because Diversatek Healthcare has no control over the conditions of use, patient selections, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, Diversatek Healthcare does not warrant either a good effect, or against an ill effect following its use. Any misuse or mishandling of the InnerVision Transillumination Fiber Optic Cable will render the foregoing limited warranty void and of no effect. If Diversatek Healthcare determines that the InnerVision Transillumination Fiber Optic Cable did not comply with the foregoing limited warranty at the time of shipment, Diversatek Healthcare will replace the nonconforming InnerVision Transillumination Fiber Optic Cable free of charge. The remedy of replacement is the sole and exclusive remedy for breach of the foregoing limited warranty. THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT.

THE FOREGOING WARRANTY WILL BE ACTIVATED AND BECOME EFFECTIVE ONLY UPON THE PURCHASER’S COMPLETION AND SUBMISSION, AND DIVERSATEK HEALTHCARE RECEIPT OF, THE APPLICABLE WARRANTY REGISTRATION CARD INCLUDED WITH THE InnerVision Transillumination Fiber Optic Cable, WITHIN THIRTY (30) CALENDAR DAYS AFTER THE PURCHASER’S RECEIPT OF THE InnerVision Transillumination Fiber Optic Cable.

**REORDERING INFORMATION**

Part No.	InnerVision Transillumination System Components
1450-40	InnerVision Fiber Optic Cable 40 Fr. Circon-ACMI/Olympus Connector (Clean, Non-sterile)
1450-50	InnerVision Fiber Optic Cable 50 Fr. Circon-ACMI/Olympus Connector (Clean, Non-sterile)
1450-56	InnerVision Fiber Optic Cable 56 Fr. Circon-ACMI/Olympus Connector (Clean, Non-sterile)
1450-02	Wolf/Dyonics Adapter (used with the InnerVision Cable with a Circon-ACMI/Olympus Connector) (Clean, Non-sterile)
1455-40	InnerVision Fiber Optic Cable 40 Fr. Storz Connector (Clean, Non-sterile)
1455-50	InnerVision Fiber Optic Cable 50 Fr. Storz Connector (Clean, Non-sterile)
1455-56	InnerVision Fiber Optic Cable 56 Fr. Storz Connector (Clean, Non-sterile)
1465-40	InnerVision Detachable Tips 40 Fr., Box of 5 Sterile
1465-50	InnerVision Detachable Tips 50 Fr., Box of 5 Sterile
1465-56	InnerVision Detachable Tips 56 Fr., Box of 5 Sterile

InnerVision® is a registered trademark of Diversatek Healthcare, Inc.

Diversatek™ is a trademark of Diversatek, Inc.

STERRAD® is a registered trademark of Advanced Sterilization Products, a division of Johnson & Johnson Medical, Inc

Detach on this line and return.

---

**Diversatek Healthcare Inc.  
InnerVision® Transillumination System  
Registration/Warranty Return Card**

THE FOREGOING WARRANTY WILL BE ACTIVATED AND BECOME EFFECTIVE ONLY UPON THE PURCHASER'S COMPLETION AND SUBMISSION, AND DIVERSATEK HEALTHCARE RECEIPT OF, THIS WARRANTY REGISTRATION CARD WITHIN THIRTY (30) CALENDAR DAYS AFTER THE PURCHASER'S RECEIPT OF THE InnerVision Transillumination Fiber Optic Cable.

Please print.

Lot Number: \_\_\_\_\_

Date: \_\_\_\_\_

Institution: \_\_\_\_\_

Dept.: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_

Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Contact Persons:

Nurse Manager: \_\_\_\_\_

Tel: \_\_\_\_\_

Other: \_\_\_\_\_

Title: \_\_\_\_\_

Tel: \_\_\_\_\_

Used by:

General Surgeon

GI Surgeon

Other: \_\_\_\_\_

Primary procedure for which this InnerVision Fiber Optic Cable is used:

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