



















**Clean**  
Approved Chemistries: Enzol® (Advanced Sterilization Products), Endozime® AW Plus (Ruhof Corporation), EmPower™ (Metrex Research), Revital-Ox™ (Steris)

1. Rinse	2. Soak	3. Wipe	4. Rinse
  <p>1 min</p> <p>Immediately after use:</p> <ol style="list-style-type: none"> <li>Install handle cap and hand-tighten securely.</li> <li>Rinse entire probe with potable water for 1 minute (Max 50°C).</li> </ol>	  <p>5 min – 30 min</p> <p>Immerse entire probe in cleaning solution and soak for 5 to 30 minutes.</p>	 <p>Detergent wipe</p> <p>Toward the end of the soak time, gently wipe sensors with gauze wet with detergent.</p>	  <p>1 min</p> <ol style="list-style-type: none"> <li>Rinse probe with potable water for 1 minute.</li> <li>Carefully dry with soft gauze or lint-free cloth.</li> <li>Visually inspect probe for any damage.</li> </ol>

**Disinfect**  
Approved Chemistries: Cidex® OPA (Advanced Sterilization Products), Cidex® 2.4% Glutaraldehyde (Advanced Sterilization Products), Revital-Ox® RESERT® (Steris)

5. Soak			6. Rinse	7. Wipe	8. Store
<p><u>Cidex® OPA</u> (non-glutaraldehyde)</p>   <p>12-30 min</p> <p>Immerse entire probe in disinfecting solution and soak for 12 to 30 minutes.</p>	<p><u>Cidex® 2.4%</u> Glutaraldehyde</p>   <p>45-60 min</p> <p>Immerse entire probe in disinfecting solution and soak for 45 to 60 minutes.</p>	<p><u>Revital-Ox® RESERT®</u></p>   <p>8-30 min</p> <p>Immerse entire probe in disinfecting solution and soak for 8 to 30 minutes.</p>	  <p>1 min</p> <p>Rinse entire probe with potable water or distilled water for a minimum of 1 minute.</p>	 <p>Dry wipe</p> <ol style="list-style-type: none"> <li>Carefully dry the catheter with a soft gauze or lint-free cloth.</li> <li>Air dry for at least 5 minutes before storing.</li> </ol>	 <p>Store the probe in a manner that keeps it dry and safe from physical damage per your institution's practices for high-level disinfected devices.</p>

 For detailed instructions and additional information refer to HRIM-195 High Resolution Impedance Manometry Probe Installation and User's Guide.

HRIM-196 Rev D, ECO 223-004

**Diversatek Healthcare**  
Technical Research & Training Center  
9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

P 800.558.6408 or 303.470.7020  
technicalsupport@diversatekhc.com  
clinicalsupport@diversatekhc.com