M-Flex[®] Blue Silicone Maloney & Hurst Esophageal Bougies

Instructions for Use and Cleaning

Read carefully prior to use.

Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

063-0026W Rev.I, 2024/02, ECO 2317

Diversatek Healthcare

Corporate Headquarters

102 E. Keefe Ave.

Milwaukee, WI 53212 USA

P 800.558.6408 | 414.265.7620

sales@diversatekhc.com

DiversatekHealthcare.com

M-FLEX® BLUE SILICONE MALONEY & HURST BOUGIES

Instructions for Use and Cleaning
Read carefully prior to use.

DESCRIPTION**M-FLEX® BLUE SILICONE TUNGSTEN MALONEY & HURST BOUGIES**

This package contains one or more blue silicone Bougies. Maloney Bougies have a tapered distal end. Hurst Bougies have a blunt distal end.

M-Flex® Blue Silicone Bougies are available in 21 sizes from 20 French to 60 French, in 2 French increments (e.g. 20FR, 22FR, 24FR...). All sizes can be purchased individually, as a complete set (contains all 21 sizes), or as a mini-set (contains 10 sizes).

INTENDED USE

For dilatation of upper esophageal webs, lower esophageal rings, caustic strictures, peptic esophageal strictures, and temporary ease of esophageal carcinoma.

INSPECTION – PRIOR TO USE

Carefully inspect each M-Flex® Blue Silicone Bougie prior to use:

1. Inspect new Bougies for damage from shipping and handling.
2. Inspect older Bougies for wear. To do this, bend and rotate the Bougie at the distal tip, shaft and proximal end. Worn Bougies should NOT be used. Signs of wear can include: surface cracks, cuts, permanently curved distal tip, and extremely faded print. Pay careful attention to bonded areas during this inspection.
3. Make sure the three (3) year expiration date has not been exceeded.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic strictures; inability to advance the bougie through the strictured area; coagulopathy; known or suspected perforation; severe inflammation or scarring near the dilation site, recent myocardial infarction, active ulcer and severe cervical arthritis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper gastrointestinal endoscopy and esophageal dilation include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

INSTRUCTIONS FOR USE

Bougies should be washed with sterile water before using in order to clean all particulates from the surface. Lubrication can be performed by wetting the surface with sterile water or lubricating with a standard medical grade water based lubricant. It is recommended that this procedure be performed under fluoroscopy. M-Flex® Blue Silicone Bougies are completely radiopaque.

1. Perform screening endoscopy to identify strictured area.
2. Generously lubricate bougie and advance it to strictured area.
Note: M-Flex® Blue Silicone Bougies have two (2) series* of numbers printed on each bougie on opposing sides. One series measures the distance to the distal tip, identified as **FROM TIP (American)**. The other series measures the distance to the stricture, identified as **FROM MAX O.D. (European)**. These printed distance markings can be used as a guide to avoid over-extending into the gastric cavity.
3. Proceed with esophageal dilation. Determine the bougie to be used by the size of the strictures, starting with a smaller French sized bougie and gradually increasing the French size. Hold the bougie at the blunt, proximal end.

4. Upon completion of esophageal dilation, remove bougie from patient.
Note: Successive dilations can be accomplished by repeating steps 2-4.
5. Clean and prepare bougie for future use following cleaning instructions below.

*Note: The Hurst style Bougie has only one series of printed numbers, FROM TIP (American).

CLEANING INSTRUCTIONS

CAUTION:

1. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, excessive wear, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.
2. **DO NOT USE** Betadine to clean the bougies as it may cause discoloration of the material.
3. **DO NOT USE** Ethylene Oxide (EO) gas sterilization.
4. **DO NOT USE** alcohol for cleaning or disinfection.
5. Appropriate gloves, gowns, and other protective wear should be worn when cleaning and sterilizing soiled or used bougies.
6. It is the responsibility of the user to qualify any deviations from the recommended methods of cleaning.

Bougies must be thoroughly cleaned prior to use.

PRE-CLEANING:

1. Immediately after use, soak bougie in a solution of warm water and enzymatic cleaning solution. **Note:** Refer to the enzymatic cleaning solution instructions for recommended duration of cleaning.
2. Remove gross contaminants on the external surfaces of the bougie by hand, using a soft bristle brush, sponge, or wipe, with warm water and enzymatic solution.
3. Thoroughly rinse the bougie with clean water.
4. Pat dry with clean, non-contaminating wipes.
5. Proceed with high-level disinfection / sterilization of bougie.

HIGH-LEVEL DISINFECTION / STERILIZATION:

The following disinfection solutions can be used with the M-Flex® Blue Silicone Bougies:

- Glutaraldehyde
- Ortho-phthalaldehyde (OPA)
- Peracetic acid
- Hydrogen peroxide

Please refer to the manufacturer's instructions for complete information and precautions for disinfection solutions with any of the above active ingredients.

The bougies have additionally been validated with the following specific brands:

- Cidex® OPA High Level Disinfectant
- Steris Resert™ XL HLD High Level Disinfectant
- Steris® V-PRO Low Temperature Sterilization System
- Rapicide PA High Level Disinfectant
- Medivators Advantage Plus System
- Medivators DSD series of automated endoscope reprocessors
- STERRAD NX® (Standard) System. The STERRAD NX® sterilization validation is based on a standard 28 minute cycle. The maximum load limit for the STERRAD NX® System is 4850 grams (10.7 lbs) per shelf. Do not use sterilization pouches/peel pouches when performing the STERRAD NX® sterilization process. Doing so may result in damage to the device. It is recommended to use polypropylene sterile wrapped trays.
- Sterrad® 100S (Short) System

Please refer to the manufacturer's instructions for complete information and precautions.

Note: Proper rinsing with clean water and complete drying of the bougie is necessary with any disinfection solution.

The M-Flex® Blue Silicone Bougies can be steam sterilized by autoclave using a gravity-type steam sterilization cycle of 132°C for 15 minutes or a prevacuum steam sterilization cycle of 132°C for 4 minutes. The Bougies should be wrapped for the steam sterilization cycle. Minimum dry time is 20 minutes. After steam sterilization, allow the Bougies to completely cool to room temperature (25°C). Do not immerse in cold liquids. Dry thoroughly.

STORAGE

The M-Flex® Blue Silicone Bougies should be stored out of direct light, and at room temperature. M-Flex® Blue Silicone Bougies should be stored flat whenever possible.

DO NOT hang the M-Flex® Blue Silicone Bougies during storage.

The Bougie Storage Case provides a protected environment.

Note: Any serious incidents that occur in relation to the M-Flex® Blue Silicone Bougies should be reported to Diversatek Healthcare and the competent authority of the Member State in which the user and/or patient is established.

DESCRIPTION OF SYMBOLS USED ON LABELS

Manufacturer

Use-by date
(YYYY-MM-DD)

Lot number

Part number

Non-sterile

Consult instructions for use


Caution, consult accompanying documents

Not made with natural rubber latex

Medical device

Authorized Representative in the European
Community

CE marking

Unique Device Identification

Diversatek™ Healthcare is a trademark of Diversatek Inc.
M-Flex® is a registered trademark of Diversatek Healthcare.

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Diversatek™ Healthcare is een handelsmerk van Diversatek Inc.