

Capture™ Polyp Trap



The smart solution to capturing,
collecting and transporting samples.

At Diversatek Healthcare, we're committed to making products that help simplify complicated procedures. By combining top-of-the-line elements with a user-friendly design, the new Capture™ Polyp Trap does just that.

Capture™ Polyp Trap

Simplify the Process.

The Capture™ Polyp Trap offers a range of features that makes the process of capturing, collecting and transporting specimens easier. Use one product to collect more polyps with less effort.

Effortless Capture.

The transparent, fixed chambers with built-in strainers feature large, easy-to-read numbers and unobstructed visibility. When suctioning, simply align the integral lid tab with the numbered chamber to capture the specimen.

Convenient Transport.

Connect the PVC tube to the suction connector to seal the Capture™ Polyp Trap for safe transport and handling post procedure.

Capture™ Polyp Trap

Single-use, non-sterile

Part No.	Description	Qty/Case
1202-01	Polyp Trap—4 Chamber	25



Four chambers with built-in strainers collect and separate specimens.



Easy-to-read numbers label each chamber.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

Learn more from
your personal product
representative at
800-558-6408