

# CareGuard® Bite Blocks



**The bold choice for upper GI procedures.**

Why compromise sensitive, expensive equipment when a simple solution from Diversatek Healthcare can change everything? Our line of bite blocks takes a colorful twist on caring for your endoscopes and your patients.

# CareGuard® Bite Blocks



## The bold choice for upper GI procedures.

You know Diversatek Healthcare for our innovative, reliable and cost-effective solutions. Now, there are even more reasons to smile. Our line of bite blocks takes a colorful twist on caring for your endoscopes and your patients. Select from 48 FR and 60 FR sizes, as well as a 60 FR XL model that features more substantial lip guards. The durable construction of CareGuard® Bite Blocks provide the protection you need.

### Features & Benefits

- 60 FR model easily accommodates a 60 FR Diversatek Healthcare Esophageal Bougie or Wire-Guided Dilator
- Dental retention rim maintains the position of the bite block
- Individually wrapped, color-coded
- 60 FR XL lateral ports allow for easy positioning or additional suctioning equipment

## CareGuard® Bite Block

Single-use, non-sterile

Part No.	Size	Description	Strap Option	Qty/Case
1196-48	48 FR	With dental rim	Blue Strap	100
1196-60	60 FR	With dental rim and lateral ports	Blue Strap	100
1196-60XL	60 FR XL	With dental rim and large flange	Blue Strap	100
1197-48	48 FR	With dental rim	N/A	100
1197-60	60 FR	With dental rim and lateral ports	N/A	100
1197-60XL	60 FR XL	With dental rim and large flange	N/A	100
1198-48	48 FR	With dental rim	Soft Strap	100
1198-60	60 FR	With dental rim and lateral ports	Soft Strap	100
1198-60XL	60 FR XL	With dental rim and large flange	Soft Strap	100



Adjustable blue strap (top), adjustable soft strap (bottom).



48FR



60FR



60FR XL

[DiversatekHealthcare.com](http://DiversatekHealthcare.com)

P 800.558.6408 or 414.265.7620

**Corporate Headquarters**  
102 East Keefe Avenue  
Milwaukee, WI 53212 USA

orders@diversatekhc.com  
sales@diversatekhc.com

**Technical Research & Training Center**  
9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

Learn more from  
your personal product  
representative at  
**800-558-6408**