


CleanFreak® Endoscope Transport Pad



Simplify transport from procedure to reprocessing.

The CleanFreak® Endoscope Transport Pad is a convenient post-procedure solution for handling endoscopes and other accessories.

CleanFreak[®] Endoscope Transport Pad

Convenient Design.

The CleanFreak[®] Endoscope Transport Pad quickly unfolds to 40" x 29" with a flexible drape design suitable for use on large and small surfaces.

Safe, Smart Transport.

After the procedure, cinching the strings on both ends of the pad will convert it into a leak-proof bag for transporting equipment for reprocessing. The plastic liner features absorbent padding to contain fluids, and is sized generously to ensure the endoscope is not coiled too tightly during transport.

CleanFreak[®] Endoscope Transport Pad Single-use, non-sterile

Part No.	Size	Qty/Case
1201-90	40" x 29"	25

Supports SGNA Guidelines* for Infection Control

- Transport the soiled endoscope to the reprocessing area in a closed container
- Containers, sinks, and basin should be large enough that the endoscope will not be damaged by being coiled too tightly

*Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes, SGNA © 2016



Edge-to-edge absorbent padding.



Pull cords around dirty endoscope to transport.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

Learn more from
your personal product
representative at
800-558-6408