


Lasso™ Rotatable Polypectomy Snare



True 360° rotation simplifies control and precision.

Lasso™ Rotatable Polypectomy Snares from Diversatek Healthcare are designed to provide ease of use and accurate placement while accommodating both hot and cold polypectomy.

Lasso™ Rotatable Polypectomy Snare

Our Smart Handle means responsive, precise control.

- Rotate the positioning barrel left or right—a full 360° with 1:1 rotation—without removing your other hand from the finger loops
- Shortened throw for controlled actuation and ease of use


Versatility in design.

- Offered in both Hot/Cold and dedicated Cold designs to meet all of your needs
- Balanced wire designs maximize tissue contact and provide a tactile closure sensation, regardless of hot or cold application


Lasso™ Rotatable Polypectomy Snares

Single-use, Sterile

Part No.	Description	Indication	Length (mm)	Minimum Working Channel (mm)	Unit
1180-01	Oval 10mm	Hot/Cold	230	2.8	10/cs
1180-02	Oval 15mm	Hot/Cold	230	2.8	10/cs
1180-03	Oval 25mm	Hot/Cold	230	2.8	10/cs
1180-04	Oval 32mm	Hot/Cold	230	2.8	10/cs
1180-10	Hex 15mm	Hot/Cold	230	2.8	10/cs
1180-11	3-in-1 Hex 10/15/25	Hot/Cold	230	2.8	10/cs
1180-12	Oval 10mm	Cold	230	2.8	10/cs
1180-13	Oval 15mm	Cold	230	2.8	10/cs
1180-14	Hex 10mm	Cold	230	2.8	10/cs
1180-15	Hex 15mm	Cold	230	2.8	10/cs
4525	Universal Active Cord	—	—	—	1 ea



Easy and responsive 360° and 1:1 rotation without the need to remove your hand.



3-in-1 Hex Snare (1180-11) is adjustable to 3 distinct diameters to accommodate multiple polyp sizes, reducing the total number of snares needed in a procedure.



Our Lasso™ Rotatable Cold Snares provide exceptional tactile grasping of polyps for precise capture and smooth cutting, even in the retroflex position.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408