

# BullDog® Disposable Endoscope Valves and Accessories



**Simplify endoscope use and clean up  
with disposable valves and accessories.**

Cleaning an endoscope properly after a procedure is a critical part of the process. Cleaning simple components can be time consuming and difficult, and also presents the risk of cross-contamination. Reduce the risk with a comprehensive portfolio of economical, single-use solutions from Diversatek Healthcare.

# Bulldog® Disposable Endoscope Valves and Accessories

## Single-use components save time and reduce risk.

Reusable endoscope valves are difficult to clean effectively and may increase procedure time and risk. Bulldog® Disposable Valves emulate the feel and action of traditional OEM valves.

## Engineered for a perfect fit.

Bulldog® Biopsy Valves accommodate larger instruments by simply lifting the cap to reveal a larger entry port. The disposable Air/Water Channel Cleaning Adapter eliminates the need to reprocess and track to the originating endoscope. And for convenience, color coding by brand ensures you'll select the right valve for Olympus®, Pentax® or Fujinon® endoscopes.



Air/Water Valves and Suction Valves feature a low-profile, responsive design for improved comfort. Compatible with Olympus, Fujinon and Pentax endoscopes.

## Bulldog® Disposable Endoscope Valves and Accessories

Part No.	Description	Unit
1195-22	Auxiliary Water Jet Connector, Olympus Compatible*, Sterile	100/cs
1203-00	Biopsy Valve, Olympus/Fujinon Compatible, Non-Sterile	100/cs
1203-05	Biopsy Valve, Olympus/Fujinon Compatible, Non-Sterile, Bulk Packed	200/cs
1203-20	Biopsy Valve, Pentax Compatible, Non-Sterile	100/cs
1203-50	Irrigating Adapter, Non-Sterile	50/cs
1203-60	Air/Water Channel Cleaning Adapter, Olympus Compatible, Non-Sterile	100/cs
1203-70	2-Piece Kit: Air/Water Valve, Suction Valve, Olympus Compatible, Sterile	50/cs
1203-71	3-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Olympus Compatible, Sterile	50/cs
1203-72	4-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Auxiliary Water Jet Connector, Olympus Compatible*, Sterile	50/cs
1203-73	5-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Auxiliary Water Jet Connector, Air/Water Cleaning Adapter, Olympus Compatible*, Sterile	40/cs
1203-81	3-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Pentax Compatible, Sterile	50/cs
1203-82	4-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Auxiliary Water Jet Connector, Pentax Compatible, Sterile	50/cs
1203-91	3-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Fujinon Compatible, Sterile	50/cs
1203-92	4-Piece Kit: Air/Water Valve, Suction Valve, Biopsy Valve, Auxiliary Water Jet Connector, Fujinon Compatible, Sterile	50/cs

\* Not compatible with Olympus® 140 Series GI Endoscopes.  
For Olympus® 140/ 160/ 180/ 190/ 240/ 260/290 Series GI Endoscopes. Olympus is a registered trademark of Olympus Optical Co., Ltd.

For Fujifilm® 700 Series Endoscopes. Fujinon is a registered trademark of Fujinon Inc.

For Pentax® 90/110/110c Series Endoscopes. Pentax is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha.



Air/Water Cleaning Adapter eliminates the need to reprocess and track. Clearly labeled for "pre-cleaning" process.



Multiple set configurations address the specific requirements of your department.

**DiversatekHealthcare.com**

P 800.558.6408 or 414.265.7620

**Corporate Headquarters**  
102 East Keefe Avenue  
Milwaukee, WI 53212 USA

orders@diversatekhc.com  
sales@diversatekhc.com

**Technical Research & Training Center**  
9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com  
clinicalsupport@diversatekhc.com  
clinicaleducation@diversatekhc.com

Learn more from  
your personal product  
representative at  
**800-558-6408**