


Capture™ Polyp Trap



Capture with Confidence.

Specimen retrieval solutions with clear visualization and uninterrupted suction make it easy to capture multiple polyps


Capture™ Polyp Trap



Capture™ 2-Tray Polyp Trap

The 2-Tray Polyp Trap features an audible and tactile, color-coded *Click-Lock* that ensures the securement of the tray.

Color-Coded Trays
track and keep multiple specimens safe, separate and secure




Tactile Click-Lock
confirms secure tray placement

Polyp Picks (2)
enable contact-free specimen retrieval

Capture™ 4-Chamber Polyp Trap

The four (4) fixed chambers with built-in strainers feature identification numbers and unobstructed visibility.

Rotatable, Easy-Lift Lid
for smooth transition between chambers



Labeled Fixed Chambers
collect four separate specimens

Innovative Strainer Design
efficiently filters specimen from fluid

Capture™ Polyp Trap

Single-use, non-sterile

Part No.	Description	Specimen Chambers/Trays	Qty/Case
1202-01	Capture 4-Chamber Polyp Trap	4	25
1202-02	Capture 2-Tray Polyp Trap	2	20

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

Learn more from
your personal product
representative at
800-558-6408