

**Zvu® Manometry Probe Reprocessing Certification**

All probes being returned to Diversatek Healthcare must be cleaned/disinfected before they will be evaluated or repaired. Any probe that is not cleaned and disinfected prior to shipment to Diversatek Healthcare will be assessed a disinfection fee.

Follow the reprocessing instructions as outlined in the IFU to clean and disinfect the Zvu Probe. If the probe is in a condition where disinfection is not possible, place the probe in a sealed bag with a biohazard label clearly affixed to the bag.

Complete this form and send with the returned probe.

<b>Probe Model Number:</b>	
<b>Probe Serial Number:</b>	

**This probe was reprocessed using:**

<u>Cleaning Solutions</u>	<u>High Level Disinfection Solutions</u>
Validated Solutions:	Validated Solutions:
<input type="checkbox"/> Revital-Ox™ Enzymatic Detergent	<input type="checkbox"/> Revital-Ox® RESERT®
<input type="checkbox"/> Enzol®	<input type="checkbox"/> Cidex® OPA
<input type="checkbox"/> Empower™	<input type="checkbox"/> Cidex® (2.4% Glutaraldehyde)
<input type="checkbox"/> Endozime® AW Plus	<i>Part# HRIM-201 Only</i>
<input type="checkbox"/> Other (list solution/method used):	<input type="checkbox"/> Other (list solution/method used):

<b>Name:</b>	
<b>Department:</b>	
<b>Organization:</b>	
<b>Phone Number:</b>	
<b>Email Address:</b>	
<b>Signature:</b>	
<b>Date:</b>	

**By signing this form, the signatory is certifying that the probe has been properly cleaned and disinfected per the Diversatek Healthcare Zvu® High Resolution Probe Reprocessing Instructions.**