

Gastroenterology Product Catalog: Devices



P 800.558.6408 or 414.265.7620

Corporate Headquarters

102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center

9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408

Table of Contents: Devices

Esophageal Dilators

M-Flex® Blue Silicone Bougies	1
WeightRight® Tungsten-Filled Bougies	1
SafeGuide® Single Disposable Over the Guidewire Esophageal Dilators	2
SafeGuide® Over the Guidewire Esophageal Dilators	2
SafeGuide® Accessories	2
Viper® 3-Stage Balloon Dilators	2
Esophageal Dilator Storage	3
InnerVision® Transillumination System	6

Hemostasis

RePlay® Hemostasis Clips	3
--------------------------------	---

Polypectomy & Tissue Acquisition

GIJaw™ Biopsy Forceps	4
Lasso™ Polypectomy Snares	4

Foreign Body Retrieval

M-Catch® Retrieval Nets	4
EndoGrip™ Grasping Forceps	4

Endoscope Care & Accessories

PureFlo® Endoscope Valves & Accessories	5
BullDog® Endoscope Valves & Accessories	5
Capture™ Polyp Trap	5
CareGuard® Bite Blocks	5

CleanFreak® Endoscope Cleaning Products

Brushes	6
Sponges & Pre-Cleaning Kit	6
Transport Pad	6
Basic Endoscopy Kits	6

Esophageal Dilators

M-Flex® Blue Silicone Bougies

M-Flex® Silicone Bougies: Maloney
#1206-00 thru 1206-60



M-Flex® Silicone Bougies: Hurst
#1207-00 thru 1207-60

M-Flex® Bougie Sets

Part No.	Description	Contents
1206-00	Maloney (tapered) Full Set	21-piece set, 20–60 FR
1206-10	Maloney (tapered) Mini Set	10-piece set, 36–54 FR
1207-00	Hurst (blunt) Full Set	21-piece set, 20–60 FR
1207-10	Hurst (blunt) Mini Set	10-piece set, 36–54 FR
1210-01	Bougie Storage Case	

M-Flex® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1206-20	1207-20	20	6.6
1206-22	1207-22	22	7.4
1206-24	1207-24	24	8
1206-26	1207-26	26	8.6
1206-28	1207-28	28	9.3
1206-30	1207-30	30	10
1206-32	1207-32	32	10.7
1206-34	1207-34	34	11.2
1206-36	1207-36	36	12
1206-38	1207-38	38	12.5

M-Flex® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1206-40	1207-40	40	13
1206-42	1207-42	42	14
1206-44	1207-44	44	14.5
1206-46	1207-46	46	15.2
1206-48	1207-48	48	16
1206-50	1207-50	50	16.5
1206-52	1207-52	52	17.2
1206-54	1207-54	54	18
1206-56	1207-56	56	18.5
1206-58	1207-58	58	19.2
1206-60	1207-60	60	20

WeightRight® Tungsten-Filled Bougies

WeightRight® Tungsten-Filled Bougies: Maloney
#1208-00 thru 1208-60



WeightRight® Tungsten-Filled Bougies: Hurst
#1209-00 thru 1209-60

WeightRight® Bougie Sets

Part No.	Description	Contents
1208-00	Maloney (tapered) Full Set	23-piece set, 16–60 FR
1208-10	Maloney (tapered) Mini Set	10-piece set, 36–54 FR
1209-00	Hurst (blunt) Full Set	23-piece set, 16–60 FR
1209-10	Hurst (blunt) Mini Set	10-piece set, 36–54 FR
1210-01	Bougie Storage Case	

WeightRight® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1208-16	1209-16	16	5.3
1208-18	1209-18	18	6
1208-20	1209-20	20	6.6
1208-22	1209-22	22	7.4
1208-24	1209-24	24	8
1208-26	1209-26	26	8.6
1208-28	1209-28	28	9.3
1208-30	1209-30	30	10
1208-32	1209-32	32	10.7
1208-34	1209-34	34	11.2
1208-36	1209-36	36	12
1208-38	1209-38	38	12.5

WeightRight® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1208-40	1209-40	40	13
1208-42	1209-42	42	14
1208-44	1209-44	44	14.5
1208-46	1209-46	46	15.2
1208-48	1209-48	48	16
1208-50	1209-50	50	16.5
1208-52	1209-52	52	17.2
1208-54	1209-54	54	18
1208-56	1209-56	56	18.5
1208-58	1209-58	58	19.2
1208-60	1209-60	60	20

Esophageal Dilators

SafeGuide® Single Disposable Over the Guidewire Esophageal Dilators

SafeGuide® Single Dilators
#1215-15 thru 1215-60



SafeGuide® Single Dilators


Part No.	French Size (FR)	Diameter (mm)
1215-15	15	5
1215-18	18	6
1215-21	21	7
1215-24	24	8
1215-27	27	9
1215-30	30	10
1215-33	33	11
1215-36	36	12

SafeGuide® Single Dilators

Part No.	French Size (FR)	Diameter (mm)
1215-39	39	13
1215-42	42	14
1215-45	45	15
1215-48	48	16
1215-51	51	17
1215-54	54	18
1215-57	57	19
1215-60	60	20

SafeGuide® Over the Guidewire Esophageal Dilators

SafeGuide® Reusable Dilators
#1214-15 thru 1214-60



SafeGuide® Reusable Dilators

Part No.	French Size (FR)	Diameter (mm)
1214-15	15	5
1214-18	18	6
1214-21	21	7
1214-24	24	8
1214-27	27	9
1214-30	30	10
1214-33	33	11
1214-36	36	12

SafeGuide® Reusable Dilators

Part No.	French Size (FR)	Diameter (mm)
1214-39	39	13
1214-42	42	14
1214-45	45	15
1214-48	48	16
1214-51	51	17
1214-54	54	18
1214-57	57	19
1214-60	60	20

SafeGuide® Dilatation Systems

Part No.	Description	Contents
1214-00	Full Set	16 dilators, 15–60 FR Includes 1 case of disposable marked spring tip guidewires, 1 case of disposable cleaning brushes, 25 disposable cleaning adapters

SafeGuide® Accessories

SafeGuide® Accessories

Part No.	Description	Qty/Case
1214-01	SafeGuide Storage Case	1
1214-02	SafeGuide Reusable Spring Tip Guidewire	2
1214-02D	SafeGuide Disposable Spring Tip Guidewire	6
1214-03D	SafeGuide Disposable Cleaning Brush	3
1214-03R	SafeGuide Reusable Cleaning Brush	3
1214-04	SafeGuide Disposable Cleaning Adapters	25



SafeGuide® Disposable
Spring Tip Guidewire
#1214-02D

Viper® 3-Stage Balloon Dilators



Viper® 3-Stage Fixed Wire Balloon Dilator
#1204-08 thru 1204-20

VIPER®



Viper® 3-Stage Wire Guided Balloon Dilator
#1205-08 thru 1205-20

Viper® 3-Stage Fixed Wire Balloon Dilators

Balloon Length 8 cm, Working Length 180 cm

Part No.	Diameter (mm)	French Size (FR)	ATM Pressure
1204-08	6-7-8	18-21-24	4-6-8
1204-10	8-9-10	24-27-30	4-6-8
1204-12	10-11-12	30-33-36	2-4-6
1204-15	12-13.5-15	36-40.5-45	2-4-6
1204-18	15-16.5-18	45-49.5-54	2-4-6
1204-20	18-19-20	54-57-60	2-4-6

Viper® 3-Stage Wire Guided Balloon Dilators

Balloon Length 5.5 cm, Working Length 240 cm

Part No.	Diameter (mm)	French Size (FR)	ATM Pressure
1205-08	6-7-8	18-21-24	4-6-8
1205-10	8-9-10	24-27-30	4-6-8
1205-12	10-11-12	30-33-36	2-4-6
1205-15	12-13.5-15	36-40.5-45	2-4-6
1205-18	15-16.5-18	45-49.5-54	2-4-6
1205-20	18-19-20	54-57-60	2-4-6

Viper® Dilatation Syringe/Gauge

Part No.	Fluid Capacity	Qty/Case
1204-01	60 cc/mL	5/cs


Esophageal Dilators

Esophageal Dilator Storage

Esophageal Dilator Storage

Includes fan with HEPA filter and 405 nm LED disinfection lighting system

Part No.	Description	Dimensions
1214-70	Esophageal Dilator Wall-mount Storage Cabinet	19"W x 4"D x 38"H
1214-71	Esophageal Dilator Horizontal Storage Cart	41"W x 23"D x 36"H
1214-72	Esophageal Dilator Vertical Storage Cart	24"W x 44"D x 24"H



Wall-Mount Storage Cabinet
#1214-70



Horizontal Storage Cart
#1214-71



Vertical Storage Cart
#1214-72

Replacement Parts & Accessories

Part No.	Description	Qty/Case
1214-73	HEPA Filter	1
1214-74	Storage Cabinet Cleaning Brush	10
1214-75	Soft Grip Gasket	1
1214-76	Fan	1
1214-77	Power Supply	1
1214-78	LED Light Strip, Single Side	1
1214-79	Swing Bar	1
1214-80	Dilator Storage Tray. Stores up to 12 bougies or dilators.	1


Hemostasis

Hemostasis

RePlay® Rotatable Hemostasis Clips

Single-use, Sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1170-01	11 mm width	230	2.8	5
1170-02	16 mm width	230	2.8	5
1170-03	18 mm width	230	2.8	5



RePlay® Rotatable Hemostasis Clip
#1170-01 thru 1170-03

Polypectomy & Tissue Acquisition

Biopsy Forceps

GIJaw™ Biopsy Forceps

Single-use, Sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1190-01	Serrated with Needle	230	2.8	10
1190-02	Serrated without Needle	230	2.8	10
1190-04	Jumbo, Serrated with Needle	230	3.2	10
1190-10	Smooth without Needle	160	2	10




Polypectomy Snares

Lasso™ Rotatable Polypectomy Snares

Single-use, Sterile

Part No.	Description	Indication	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1180-01	Oval 10 mm	Hot/Cold	230	2.8	10
1180-02	Oval 15 mm	Hot/Cold	230	2.8	10
1180-03	Oval 25 mm	Hot/Cold	230	2.8	10
1180-04	Oval 32 mm	Hot/Cold	230	2.8	10
1180-10	Hex 15 mm	Hot/Cold	230	2.8	10
1180-11	3-in-1 Hex 10/15/25 mm	Hot/Cold	230	2.8	10
1180-12	Oval 10 mm	Cold	230	2.8	10
1180-13	Oval 15 mm	Cold	230	2.8	10
1180-14	Hex 10 mm	Cold	230	2.8	10
1180-15	Hex 15 mm	Cold	230	2.8	10
4525	Universal Active Cord	—	—	—	1



Foreign Body Retrieval

Retrieval Nets

M-Catch® Retrieval Nets

Single-use, Non-sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1175-00	Standard Oval, 3 x 6 cm	160	2.8	5
1175-01	Standard Oval, 3 x 6 cm	230	2.8	5
1175-02	Standard Oval, 3 x 6 cm	350	2.8	5
1175-10	Polyp Oval, 3 x 6 cm	230	2.8	5




Grasping Forceps

EndoGrip™ Rotatable Grasping Forceps

Single-use, Non-sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1177-01	Rat Tooth-Alligator Combo	230	2.8	5



Endoscope Care & Accessories

Disposable Endoscope Valves & Accessories

PureFlo® Endoscope Valves & Accessories

Single-use

Part No.	Description	Piece/Kit	Compatibility	Sterility	Qty/Case
1195-21	Auxiliary Water Jet Connector*		Olympus	Sterile	100



PureFlo® Auxiliary Water Jet Connector #1195-21

BullDog® Endoscope Valves & Accessories

Single-use

Part No.	Description	Piece/Kit	Compatibility	Sterility	Qty/Case
1195-22	Auxiliary Water Jet Connector*		Olympus	Sterile	100
1203-00	Biopsy Valve		Olympus/Fujinon	Non-sterile	100
1203-05	Biopsy Valve, bulk		Olympus/Fujinon	Non-sterile	200
1203-20	Biopsy Valve		Pentax	Non-sterile	100
1203-50	Irrigating Adapter			Non-sterile	50
1203-60	Air/Water Channel Cleaning Adapter		Olympus	Non-sterile	100
1203-70	Air/Water Valve, Suction Valve	2-Piece Kit	Olympus	Sterile	50
1203-71	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Olympus	Sterile	50
1203-72	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector*, Biopsy Valve	4-Piece Kit	Olympus	Sterile	50
1203-73	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector*, Biopsy Valve, Air/Water Cleaning Adapter	5-Piece Kit	Olympus	Sterile	40
1203-81	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Pentax	Sterile	50
1203-82	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector, Biopsy Valve	4-Piece Kit	Pentax	Sterile	50
1203-91	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Fujinon	Sterile	50
1203-92	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector, Biopsy Valve	4-Piece Kit	Fujinon	Sterile	50



Olympus® 5-Piece Kit #1203-73



Pentax® 4-Piece Kit #1203-82

* Not compatible with Olympus® 140 Series GI Endoscopes.

For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes. Olympus is a registered trademark of Olympus Optical, Co., Ltd.

For Fujifilm® 700 Series Endoscopes. Fujinon is a registered trademark of Fujinon, Inc.

For Pentax® 90/110/110c Series Endoscopes. Pentax is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha.

Polyp Trap

Capture™ Polyp Trap

Single-use, Non-sterile

Part No.	Description	Qty/Case
1202-01	4 Chamber	25
1202-02	2-Tray	20



Capture® Polyp Trap: 4 Chamber #1202-01



Capture® Polyp Trap: 2-Tray #1202-02

Bite Blocks

CareGuard® Bite Blocks

Single-use, Non-sterile

Part No.	Size	Color	Style	Strap Option	Qty/Case
1196-48	48 FR	Green	With dental rim	Blue Strap	100
1196-60	60 FR	Pink	With dental rim and lateral ports	Blue Strap	100
1196-60XL	60 FR XL	Purple	With dental rim and lateral ports	Blue Strap	100
1197-60	60 FR	Pink	With dental rim and lateral ports	No Strap	100
1197-60XL	60 FR XL	Purple	With dental rim and lateral ports	No Strap	100
1198-48	48 FR	Green	With dental rim	Soft Strap	100
1198-60	60 FR	Pink	With dental rim and lateral ports	Soft Strap	100
1198-60XL	60 FR XL	Purple	With dental rim and lateral ports	Soft Strap	100



Careguard® Bite Block, 60 FR #1196-60, 1197-60, 1198-60

CleanFreak® Endoscope Cleaning Products

Cleaning Brushes

CleanFreak® Endoscope Cleaning Brushes

Single-use, Non-sterile

Part No.	Description	Bristle Diameter (mm)	Sheath Diameter (mm)	Channel Diameter Range (mm)	Length (in/cm)	Qty/Case
1200-03	Extra Long Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-04	Dual Ended Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-05	Combination Brush (Valve & Channel Brush)	5-7	1.7	2.0-6.0	90.5/230	50
1200-06	Valve Brush	10				50
1200-07	Heavy Duty Standard Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-08	Heavy Duty Dual Ended Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-09	Heavy Duty Combination Brush (Valve & Channel Brush)	5-7	1.7	2.0-6.0	90.5/230	50



CleanFreak®
Heavy Duty
Combination Brush
#1200-09

Sponges & Pre-Cleaning Kit

CleanFreak® Endoscope Cleaning Sponges

Single-use, Non-sterile

Part No.	Description	Style	Color	Qty/Case
1201-01	Endoscope Sponge	Textured	Blue	200
1201-02	Endoscope Sponge	Cylindrical	Yellow	50
1201-03	Dry Imbibed Endoscope Sponge	Cylindrical	Yellow	50

CleanFreak® Pre-Cleaning Kit

Single-use, Non-sterile

Part No.	Description	Qty/Case
1201-20	Dry Imbibed Sponge with Water Basin	24



CleanFreak®
Pre-Cleaning Kit
#1201-20

Transport Pad

CleanFreak® Transport Pad

Single-use, Non-sterile

Part No.	Description	Size	Qty/Case
1201-90	Transport Pad	40" x 29"	25



CleanFreak®
Transport Pad
#1201-90

Endoscopy Kits

Basic Endoscopy Kits

Diversatek Healthcare will customize a kit to your exact specifications. Please contact your Sales Representative for details on customizing your own kit.



Basic Endoscopy Kit

Lighted Bougies

InnerVision® Transillumination System (Lighted Bougie)

InnerVision® Fiber Optic Cable

Reusable, Non-sterile

Part No.	Size	Compatible Connector	Qty/Case
1450-40	40 Fr	Circon-ACMI/Olympus	1
1450-50	50 Fr	Circon-ACMI/Olympus	1
1450-56	56 Fr	Circon-ACMI/Olympus	1
1455-40	40 Fr	Storz	1
1455-50	50 Fr	Storz	1
1455-56	56 Fr	Storz	1

InnerVision® Fiber Optic Tips

Single-use, Sterile

Part No.	Size	Qty/Case
1465-40	40 Fr, Blunt	5
1465-50	50 Fr, Tapered	5
1465-56	56 Fr, Tapered	5

Detachable InnerVision® Tips must be used with the corresponding InnerVision® Cable




Wolf/Dyonics Adapter for Circon-ACMI/Olympus Light Source

Reusable, Non-sterile

Part No.	Qty/Case
1450-02	1

Wolf/Dyonics Adapter, Storz Connector and Circon-ACMI/Olympus Connector.

Gastroenterology Product Catalog: Testing Systems



P 800.558.6408 or 414.265.7620

Corporate Headquarters

102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center

9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408

Table of Contents: Testing Systems

Mucosal Integrity

MiVu™ Mucosal Integrity Testing System..... 1

Reflux Monitoring

ZepHr® Recorder & Accessories..... 2

ZepHr® Replacement Parts..... 2

Z/pH Probes 3

pH Probes 3

GI Motility

High Resolution Esophageal Manometry Base Design Options..... 4

High Resolution Esophageal Manometry Probes & Accessories 5

Pelvic Floor Evaluation

Anorectal Manometry Base Design Options 6

Anorectal Manometry Probes & Accessories 6–7

Computer Kits & Accessories

Replacement Computer Accessories..... 7

Product Training

Diversatek University Training 8

Mucosal Integrity

MiVu™ Mucosal Integrity Testing System for Use With MiVu™ Endo Cap

Choose a MiVu Base Design Option:

MiVu™ Base Option A: New User with No Existing Base System

For new users, who want to perform MiVu Mucosal Integrity Studies without an existing Base System

Part No.	Description
PRIZM-CART-II	PriZm® Gen-II Base Cart System. Includes PriZm® Gen-II Central Unit, 24" LCD medical grade monitor, computer, keyboard/mouse and printer, Zvu® and BioVIEW® Software.

MiVu™ Base Option B: New User with No Existing Base System. Isolation Transformer Included.

For new users, who want to perform MiVu Mucosal Integrity Studies without an existing Base System – and includes an isolation transformer

Part No.	Description
PRIZMSYS-II	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit, plus isolation transformer.

MiVu™ Base Option C: Upgrade Option for Existing inSIGHT Ultima® User

For existing inSIGHT Ultima® Base Cart System (ULTCART) or inSIGHT Ultima® Base System (ULTSYS) users who want to perform MiVu™ Mucosal Integrity Studies

Part No.	Description
PRIZMUPG-KIT	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit, Zvu® and BioVIEW® Software.


MiVu™ Endo Cap & Accessories

Part No.	Description	Qty/Case
MI-ESO-CAP-3L	MiVu™ Esophageal Endo Cap. Compatible with an endoscope diameter of 9.2 - 11.2 mm (single-use, non-sterile).	1
MI-CAB-02	MiVu™ Endo Cap Cable. Connects the MiVu™ Esophageal Endo Cap to the PriZm® Central Unit.	1

MiVu™ Mucosal Integrity Testing System is covered by one or more of the following patents:
US 11,291,382, US 9,814,408 and US 10,321,867.



PriZm® Gen-II Base Cart System
#PRIZM-CART-II



MiVu™ Esophageal Endo Cap
#MI-ESO-CAP-3L

MiVu™ Endo Cap Cable
#MI-CAB-02

Zvu® Software

GI Motility Software

Included in Base Options listed above

Part No.	Description
ZVU-3	Zvu® Software for MiVu™ Mucosal Integrity Studies, Reflux Monitoring, and Esophageal Manometry
ZVU-3-ZIP	Zvu® Software for MiVu™ Mucosal Integrity Studies, Reflux Monitoring, and Esophageal Manometry. Electronic delivery.

Reflux Monitoring (Impedance/pH)

ZepHr® Recorder & Accessories

ZepHr® Recorder

Part No.	Product	Description	Qty/Case
Z07-2000B-B	ZepHr® Z/pH Recorder	For ambulatory impedance/pH reflux monitoring. Includes carrying case, two SD cards, and two AA batteries. Limited two year warranty	1

ZepHr® Accessories

Part No.	Product	Description	Qty/Case
Z07-9000	ZepHr® System Starter Kit	Includes an SD Card Reader, pH Calibration Tubes, and Buffer Solutions. Includes Zvu® Software	1
Z07-PACK	ZepHr® Pack	ZepHr® associated accessories	1
A04-2000-S	ZepHr® Sphincter Locator	LES locator for pH or Z/pH reflux testing. Infused probes are required for use.	1
D82-3103	pH 1 Buffer Solution (500 mL)	For use with pH reflux testing system and probes	1
D82-3151	pH 4 Buffer Solution (500 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and probes	1
D82-3152	pH 7 Buffer Solution (500 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and probes	1
NG-0050	Nasal Tube Fastener	For reflux studies. Adult size	50



ZepHr®
ZepHr® Z/pH Recorder
#Z07-2000B-B



ZepHr® Replacement Products

ZepHr® Replacement Products

Part No.	Description	Qty/Case
Z07-4053	ZepHr® Battery Cover	1
Z07-4055	ZepHr® Probe Cover, Smoke	1
Z07R-8002	ZepHr® Thermoform Pouch	1
Z07-1250	SD Flash Memory Card	1
L89-5700-4	USB Multicard Reader Kit for ZepHr®	1
L89-6020-40	pH Calibration Tubes (Set of 3). Includes pH 4 and pH 7 buffer tubes, plus a water tube	1
L89-6021-40	pH Calibration Buffer Tube, 40 cm. Includes both pH 4 and pH 7 labels	1
L89-8080	pH Calibration Stand	1
L89-8081-40	pH Calibration Water Tube, 40 cm	1

High Resolution Esophageal Manometry

Choose a Manometry Base Design Option:

Manometry Base Option A: New User with No Existing Base System

For new users who want to perform High Resolution Manometry or MiVu™ Mucosal Integrity Studies, without an existing Base System

Part No.	Description
PRIZMCART-II	PriZm® Gen-II Base Cart System. Includes PriZm® Gen-II Central Unit, 24" LCD medical grade monitor, computer, keyboard/mouse and printer.
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the PriZm® Base Cart System. Includes Zvu® and BioView® Software

Manometry Base Option B: New User with No Existing Base System. Isolation Transformer Included.

For new users who want to perform High Resolution Manometry or MiVu™ Mucosal Integrity Studies, without an existing Base System – and include an isolation transformer.

Part No.	Description
PRIZMSYS-II	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit plus isolation transformer
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the PriZm® Base Cart System. Includes Zvu® and BioView® Software.

Manometry Base Option C: Upgrade Option for Existing inSIGHT Ultima® User

For existing inSIGHT Ultima® Base Cart System (ULTCART) or inSIGHT Ultima® Base System (ULTSYS) users who want to perform High Resolution Manometry and have access to MiVu™ functionality.

Part No.	Description
PRIZMUPG-KIT	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit, Zvu® and BioVIEW® Software.

Manometry Base Option D: New User with No Existing Base System

For new users who want to perform High Resolution Manometry. This option is not compatible for MiVu™ Mucosal Integrity Studies.

Part No.	Description
ULTCART	inSIGHT Ultima® Base Cart System. Includes inSIGHT Ultima® Central Unit, 24" LCD medical grade monitor, computer, keyboard/mouse and printer. Motility Cart with isolation transformer
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the inSIGHT Ultima® or PriZm® Base Cart System. Includes Zvu® and BioView® Software

Manometry Base Option E: New User with No Existing Base System

For new users who want to perform High Resolution Manometry without an existing Base System. This option is not compatible for MiVu™ Mucosal Integrity Studies.

Part No.	Description
ULTSYS	inSIGHT Ultima® Base System. Includes inSIGHT Ultima® Central Unit plus isolation transformer. • Motility Cart and computer/monitor not included
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the inSIGHT Ultima® or PriZm® Base Cart System. Includes Zvu® and BioView® Software

GI Motility Software


Included in Base Options listed above

Part No.	Description
ZVU-3	Zvu® Software for MiVu™ Mucosal Integrity Studies, Reflux Monitoring, and Esophageal Manometry
ZVU-3-ZIP	Zvu® Software for MiVu™ Mucosal Integrity Studies, Reflux Monitoring, and Esophageal Manometry. Electronic delivery.



PriZm® Gen-II
Central Unit

PriZm®



Zvu® Software
#ZVU-3

Zvu®

High Resolution Esophageal Manometry

Probe Kits for PriZm® Esophageal Manometry Systems

Motility cart, computer, monitor, and software not included

Part No.	Description	Qty/Case
HRIM-201-12-KIT *	Zvu® High Resolution Esophageal Manometry Probe Kit (34 circumferential pressure / 18 impedance channels. 12 FR. Part #HRIM-201), plus appropriate cable	1
HRIM-KIT	Unisensor HRIM2 High Resolution Impedance Manometry Probe Kit (32 circumferential pressure / 16 impedance channels. 12 FR. Part #UNI-ESO-WG1A1), plus appropriate cable.	1
HRIMSD-KIT	Unisensor HRIM SD High Resolution Impedance Manometry Probe Kit (32 directional pressure / 16 impedance channels. 8 FR. Part #UNI-ESO-WG1A2), plus appropriate cable. 12 week lead time	1

Esophageal Manometry Probes


Part No.	Description	Qty/Case
HRIM-201 *	Zvu® HRIM Esophageal Manometry Probe. 34 circumferential pressure / 18 impedance channels. 12 FR	1
UNI-ESO-WG1A1	Unisensor HRIM2 High Resolution Impedance Manometry Probe. 32 circumferential pressure / 16 impedance channels. 12 FR	1
UNI-ESO-WG1A2	Unisensor HRIM SD High Resolution Impedance Manometry Probe. 32 directional pressure / 16 impedance channels. 8 FR	1

Esophageal Manometry Accessories

Part No.	Description	Qty/Case
EFTV-12V	Viscous Solution, 80 mL containers. Used with esophageal probes utilizing impedance technology	12
H12R-5100	Air Calibration Tube Kit, 18" for inSIGHT Ultima® and PriZm®	1
H12R-5200	Water Calibration Tube Kit for inSIGHT Ultima® and PriZm®	1
H12R-5400	Auto Calibration Air Lines Kit for Zvu® Software, PriZm®, and inSIGHT Ultima®	1
H12R-7110	inSIGHT Ultima® and PriZm® Cable for Unisensor Probe	1
H12R-7610 *	Cable for Zvu® HRIM and HRaM Probes	1
ISO-MED-3	Isolation Transformer	1
S90-5120-H	Calibration Tube Universal Washer and O-Ring Replacement Kit	1



Zvu® HRIM Esophageal Manometry Probe
#HRIM-201



Zvu® Software
#ZVU-3

Zvu®

*Zvu® software version 3.0 or higher is required.

Pelvic Floor Evaluation

Anorectal Manometry

Choose a Manometry Base Design Option:

Manometry Base Option A: New User with No Existing Base System

For new users who want to perform High Resolution Manometry or MiVu™ Mucosal Integrity Studies, without an existing Base System

Part No.	Description
PRIZMCART-II	PriZm® Gen-II Base Cart System. Includes PriZm® Gen-II Central Unit, 24" LCD medical grade monitor, computer, keyboard/mouse and printer.
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the PriZm® Base Cart System. Includes Zvu® and BioView® Software



PriZm® Gen-II
Central Unit

PriZm®

Manometry Base Option B: New User with No Existing Base System. Isolation Transformer Included.

For new users who want to perform High Resolution Manometry or MiVu™ Mucosal Integrity Studies, without an existing Base System – and include an isolation transformer.

Part No.	Description
PRIZMSYS-II	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit plus isolation transformer
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the inSIGHT Ultima® or PriZm® Base Cart System. Includes Zvu® and BioView® Software

Manometry Base Option C: Upgrade Option for Existing inSIGHT® Ultima User

For existing inSIGHT Ultima® Base Cart System (ULTCART) or inSIGHT Ultima® Base System (ULTSYS) users who want to perform High Resolution Manometry and have access to MiVu™ functionality.

Part No.	Description
PRIZMUPG-KIT	PriZm® Gen-II Base System. Includes PriZm® Gen-II Central Unit, Zvu® and BioVIEW® Software.

Manometry Base Option D: New User with No Existing Base System

For new users who want to perform High Resolution Manometry. This option is not compatible for MiVu™ Mucosal Integrity Studies.

Part No.	Description
ULTCART	inSIGHT Ultima® Base Cart System. Includes inSIGHT Ultima® Central Unit, 24" LCD medical grade monitor, computer, keyboard/mouse and printer. Motility Cart with isolation transformer
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the inSIGHT Ultima® or PriZm® Base Cart System. Includes Zvu® and BioView® Software

Manometry Base Option E: New User with No Existing Base System

For new users who want to perform High Resolution Manometry without an existing Base System. This option is not compatible for MiVu Mucosal Integrity Studies.

Part No.	Description
ULTSYS	inSIGHT Ultima® Base System. Includes inSIGHT Ultima® Central Unit plus isolation transformer. Motility Cart and computer/monitor not included.
MANO-KIT	Manometry Starter Kit. For customers with an existing inSIGHT Ultima® or PriZm® System. Includes air and water calibration kits, plus connecting tubing to the inSIGHT Ultima® or PriZm® Base Cart System. Includes Zvu® and BioView® Software



Zvu® HRaM Anorectal Manometry Probe
#HRAM-200

Probe Kits for Anorectal Manometry Systems

Choose Probe Kit(s):

Part No.	Description
HRAM-200-17-KIT	Zvu® High Resolution Anorectal Manometry Probe Kit (26 unidirectional pressure / 1 directional balloon channel. 16.7 FR. Part #HRAM-200). • Plus appropriate cables, ligator, anorectal balloon procedure kit • BioVIEW® Software included • Motility Cart, computer, and monitor not included
HRAM-KIT	Unisensor High Resolution Anorectal Manometry Probe Kit (22 directional pressure / 1 directional balloon channel. 12 FR. Part #UNI-ANO-M0138). • Plus appropriate cables, ligator, anorectal balloons (10 ea), and disposable suture removal kit • Motility Cart, computer, monitor, and software not included
HRCA-KIT	Unisensor High Resolution Anorectal Manometry Probe Kit (6 circumferential pressure / 1 directional balloon channel. 12 FR. Part #UNI-ANO-704G1) • Plus appropriate cables, ligator, anorectal balloons (10 ea), and disposable suture removal kit • Motility Cart, computer, monitor, and software not included. 12 week lead time
EMG-KIT	Anal Sphincter Electromyography Kit includes EMG Adapter, cables and perianal sensor • Motility Cart, computer, monitor, and software not included

Anorectal Manometry Probes

Part No.	Description
HRAM-200	Zvu® HRaM Anorectal Manometry Probe. 26 Unidirectional pressure / 1 balloon pressure channels. 16.7 FR
UNI-ANO-704G1	Unisensor HRCA Anorectal Manometry Probe. 6 Circumferential pressure channels and 1 directional balloon channel. 12 FR. 12 week lead time
UNI-ANO-M0138	Unisensor HRaM Anorectal Manometry Probe. 22 Directional pressure channels and 1 directional balloon channel. 12 FR.

Pelvic Floor Evaluation

Anorectal Manometry

GI Motility Software

Part No.	Description	Qty/Case
S98-1925-P1	BioVIEW® Suite Software for Anorectal Studies	1
S96-1118	Data Extraction Utility Software Kit (compatible only with BioVIEW® software)	1
HRAM-100	London Protocol Update, thumb drive (compatible only with BioVIEW® software)	1
HRAM-100-ZIP	London Protocol Update, electronic delivery (compatible only with BioVIEW® software)	1

Anorectal Manometry Accessories

Part No.	Description	Qty/Case
A86-1070	Anorectal Balloon, adult, contains Latex. For anorectal compliance testing	5
A86-4100	Polyurethane Bands (disk of 20 bands), for Unisensor HRaM Probes	1
A86-4110-10	HRaM Kit. Includes 10 Anorectal Balloons, 20 Polyurethane Bands, and a Suture Removal Kit, for Unisensor HRaM Probes	1
A86-4120	HRaM Syringe/Stopcock Kit	1
A86-4200	Ligator for Unisensor HRaM Probes	1
A86-4300-10	Anorectal Balloon, 400CC TPE, for Unisensor HRaM Probes	10
A86-4500	Disposable Suture Removal Kit for Unisensor HRaM Probes	1
A86-5070	Anorectal Balloon, small. Contains Latex	5
A97-1000	StimSense Pressure Sensor. Measures rectal balloon pressure during anorectal compliance testing	1
H12R-4300	EMG Adapter Assembly. Connects to the inSIGHT Ultima® and PriZm®	1
H12R-5400	Auto Calibration Air Lines Kit for Zvu® Software, PriZm®, and inSIGHT Ultima®	1
H12R-7110	Ultima® and PriZm® Cable for Unisensor High Resolution Manometry Probe	1
H12R-7210	Ultima® and PriZm® Cable for HRCA Manometry Probe (Lemo Connector)	1
H12R-7610	Cable for Zvu® HRiM and HRaM Probes	1
HRAM-BAL-10	Anorectal Balloon for Zvu® HRaM Probes, 400CC TPE, and Balloon Retention Bands	10
HRAM-LIGATOR	Zvu® Ligator. For balloon attachment to a Zvu® HRaM Probe	1
HRAM-PROC-10	Zvu® Anorectal Balloon Attachment & Removal Kit (single-use). Includes Zvu® Anorectal Balloons, Balloon Retention Bands, Band EZ-Remove Guides, and disposable scalpel	10
ISO-MED-3	Isolation Transformer	1
SRS-4510	Perianal Rectal Sensor (single-use)	1
SR1B	Mui Scientific Anorectal Expulsion Balloon Catheter (single-use)	10

Mui Scientific Water Perfusion Catheters

Single-use. 4 Week lead time


Part No.	Description	Qty/Case
SUAMC-M5A-10	Mui Anorectal Manometric Catheter, radial configuration (14 FR, 4 channel)	10
SUAMC-M83A-10	Mui Anorectal Manometric Catheter (14 FR, 8 channel)	10
SUAMC-M8AA-10	Mui Anorectal Manometric Catheter, radial configuration (14 FR, 8 channel)	10

Computer Kits & Accessories


Replacement Computer Kits & Accessories

For inSIGHT Ultima® and PriZm® Systems

Part No.	Description	Qty/Case
COMP-KIT	Desktop PC with Windows 10. Includes a wireless keyboard and mouse • Installed with Zvu® and BioVIEW® Software	1
L89-5355	Printer USB Cable	1
L89-5500	Desktop Color Printer	1
SYS-MONITOR	Monitor, Touchscreen LCD, 24" for inSIGHT Ultima® and PriZm®	1



BioVIEW® Suite Software
#S98-1925-P1



Zvu® Anorectal Balloon Attachment & Removal Kit
(All Items Shown Above)
#HRAM-PROC-10



Zvu® Ligator
#HRAM-LIGATOR

Product Training

Product Training at Diversatek University Denver & Onsite

Product Training at Diversatek University Denver

Part No.	Description	Qty
DEN-FGID	• One Registration—3 Day Functional GI Diagnostic Training • Training on High Resolution Esophageal Manometry, Reflux Studies, and Anorectal Manometry • Location: Denver • Expires 1 year from Invoice Date (Does not include transportation/hotel costs)	1

Onsite Clinical Training

- Includes one set of master manuals
- Expires 1 year from Invoice Date

Part No.	Description	Qty
ON-SITE-ANO	Clinical Training Onsite (Two Day Session)—High Resolution Anorectal Manometry	1
ON-SITE-ESO	Clinical Training Onsite (Two Day Session)— High Resolution Esophageal Manometry	1
ON-SITE-REFLUX	Clinical Training Onsite (One Day Session)—Reflux Testing	1

Virtual Training

Part No.	Description	Qty
VIRT-FULLDAY	Full Day Virtual Training Session	1
VIRT-HALFDAY	Half Day Virtual Training Session	1

