

Gastroenterology Product Catalog: Devices



P 800.558.6408 or 414.265.7620

Corporate Headquarters

102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center

9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408

Table of Contents: Devices

Esophageal Dilators

M-Flex® Blue Silicone Bougies	1
WeightRight® Tungsten-Filled Bougies	1
SafeGuide® Single Disposable Over the Guidewire Esophageal Dilators	2
SafeGuide® Over the Guidewire Esophageal Dilators	2
SafeGuide® Accessories	2
Viper® 3-Stage Balloon Dilators	2
Esophageal Dilator Storage	3
InnerVision® Transillumination System	6

Hemostasis

RePlay® Hemostasis Clips	3
--------------------------------	---

Polypectomy & Tissue Acquisition

GIJaw™ Biopsy Forceps	4
Lasso™ Polypectomy Snares	4

Foreign Body Retrieval

M-Catch® Retrieval Nets	4
EndoGrip™ Grasping Forceps	4

Endoscope Care & Accessories

PureFlo® Endoscope Valves & Accessories	5
BullDog® Endoscope Valves & Accessories	5
Capture™ Polyp Trap	5
CareGuard® Bite Blocks	5

CleanFreak® Endoscope Cleaning Products

Brushes	6
Sponges & Pre-Cleaning Kit	6
Transport Pad	6
Basic Endoscopy Kits	6

Esophageal Dilators

M-Flex® Blue Silicone Bougies

M-Flex® Silicone Bougies: Maloney
#1206-00 thru 1206-60



M-Flex® Silicone Bougies: Hurst
#1207-00 thru 1207-60

M-Flex® Bougie Sets

Part No.	Description	Contents
1206-00	Maloney (tapered) Full Set	21-piece set, 20–60 FR
1206-10	Maloney (tapered) Mini Set	10-piece set, 36–54 FR
1207-00	Hurst (blunt) Full Set	21-piece set, 20–60 FR
1207-10	Hurst (blunt) Mini Set	10-piece set, 36–54 FR
1210-01	Bougie Storage Case	

M-Flex® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1206-20	1207-20	20	6.6
1206-22	1207-22	22	7.4
1206-24	1207-24	24	8
1206-26	1207-26	26	8.6
1206-28	1207-28	28	9.3
1206-30	1207-30	30	10
1206-32	1207-32	32	10.7
1206-34	1207-34	34	11.2
1206-36	1207-36	36	12
1206-38	1207-38	38	12.5

M-Flex® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1206-40	1207-40	40	13
1206-42	1207-42	42	14
1206-44	1207-44	44	14.5
1206-46	1207-46	46	15.2
1206-48	1207-48	48	16
1206-50	1207-50	50	16.5
1206-52	1207-52	52	17.2
1206-54	1207-54	54	18
1206-56	1207-56	56	18.5
1206-58	1207-58	58	19.2
1206-60	1207-60	60	20

WeightRight® Tungsten-Filled Bougies

WeightRight® Tungsten-Filled Bougies: Maloney
#1208-00 thru 1208-60



WeightRight® Tungsten-Filled Bougies: Hurst
#1209-00 thru 1209-60

WeightRight® Bougie Sets

Part No.	Description	Contents
1208-00	Maloney (tapered) Full Set	23-piece set, 16–60 FR
1208-10	Maloney (tapered) Mini Set	10-piece set, 36–54 FR
1209-00	Hurst (blunt) Full Set	23-piece set, 16–60 FR
1209-10	Hurst (blunt) Mini Set	10-piece set, 36–54 FR
1210-01	Bougie Storage Case	

WeightRight® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1208-16	1209-16	16	5.3
1208-18	1209-18	18	6
1208-20	1209-20	20	6.6
1208-22	1209-22	22	7.4
1208-24	1209-24	24	8
1208-26	1209-26	26	8.6
1208-28	1209-28	28	9.3
1208-30	1209-30	30	10
1208-32	1209-32	32	10.7
1208-34	1209-34	34	11.2
1208-36	1209-36	36	12
1208-38	1209-38	38	12.5

WeightRight® Maloney & Hurst Bougies

Maloney (Tapered) Part No.	Hurst (Blunt) Part No.	French Size (FR)	Diameter (mm)
1208-40	1209-40	40	13
1208-42	1209-42	42	14
1208-44	1209-44	44	14.5
1208-46	1209-46	46	15.2
1208-48	1209-48	48	16
1208-50	1209-50	50	16.5
1208-52	1209-52	52	17.2
1208-54	1209-54	54	18
1208-56	1209-56	56	18.5
1208-58	1209-58	58	19.2
1208-60	1209-60	60	20

Esophageal Dilators

SafeGuide® Single Disposable Over the Guidewire Esophageal Dilators

SafeGuide® Single Dilators
#1215-15 thru 1215-60



SafeGuide® Single Dilators


Part No.	French Size (FR)	Diameter (mm)
1215-15	15	5
1215-18	18	6
1215-21	21	7
1215-24	24	8
1215-27	27	9
1215-30	30	10
1215-33	33	11
1215-36	36	12

SafeGuide® Single Dilators

Part No.	French Size (FR)	Diameter (mm)
1215-39	39	13
1215-42	42	14
1215-45	45	15
1215-48	48	16
1215-51	51	17
1215-54	54	18
1215-57	57	19
1215-60	60	20

SafeGuide® Over the Guidewire Esophageal Dilators

SafeGuide® Reusable Dilators
#1214-15 thru 1214-60



SafeGuide® Reusable Dilators

Part No.	French Size (FR)	Diameter (mm)
1214-15	15	5
1214-18	18	6
1214-21	21	7
1214-24	24	8
1214-27	27	9
1214-30	30	10
1214-33	33	11
1214-36	36	12

SafeGuide® Reusable Dilators

Part No.	French Size (FR)	Diameter (mm)
1214-39	39	13
1214-42	42	14
1214-45	45	15
1214-48	48	16
1214-51	51	17
1214-54	54	18
1214-57	57	19
1214-60	60	20

SafeGuide® Dilatation Systems

Part No.	Description	Contents
1214-00	Full Set	16 dilators, 15–60 FR Includes 1 case of disposable marked spring tip guidewires, 1 case of disposable cleaning brushes, 25 disposable cleaning adapters

SafeGuide® Accessories

SafeGuide® Accessories

Part No.	Description	Qty/Case
1214-01	SafeGuide Storage Case	1
1214-02	SafeGuide Reusable Spring Tip Guidewire	2
1214-02D	SafeGuide Disposable Spring Tip Guidewire	6
1214-03D	SafeGuide Disposable Cleaning Brush	3
1214-03R	SafeGuide Reusable Cleaning Brush	3
1214-04	SafeGuide Disposable Cleaning Adapters	25



SafeGuide® Disposable
Spring Tip Guidewire
#1214-02D

Viper® 3-Stage Balloon Dilators



Viper® 3-Stage Fixed Wire Balloon Dilator
#1204-08 thru 1204-20

VIPER®



Viper® 3-Stage Wire Guided Balloon Dilator
#1205-08 thru 1205-20

Viper® 3-Stage Fixed Wire Balloon Dilators

Balloon Length 8 cm, Working Length 180 cm

Part No.	Diameter (mm)	French Size (FR)	ATM Pressure
1204-08	6-7-8	18-21-24	4-6-8
1204-10	8-9-10	24-27-30	4-6-8
1204-12	10-11-12	30-33-36	2-4-6
1204-15	12-13.5-15	36-40.5-45	2-4-6
1204-18	15-16.5-18	45-49.5-54	2-4-6
1204-20	18-19-20	54-57-60	2-4-6

Viper® 3-Stage Wire Guided Balloon Dilators

Balloon Length 5.5 cm, Working Length 240 cm

Part No.	Diameter (mm)	French Size (FR)	ATM Pressure
1205-08	6-7-8	18-21-24	4-6-8
1205-10	8-9-10	24-27-30	4-6-8
1205-12	10-11-12	30-33-36	2-4-6
1205-15	12-13.5-15	36-40.5-45	2-4-6
1205-18	15-16.5-18	45-49.5-54	2-4-6
1205-20	18-19-20	54-57-60	2-4-6

Viper® Dilatation Syringe/Gauge

Part No.	Fluid Capacity	Qty/Case
1204-01	60 cc/mL	5/cs


Esophageal Dilators

Esophageal Dilator Storage

Esophageal Dilator Storage

Includes fan with HEPA filter and 405 nm LED disinfection lighting system

Part No.	Description	Dimensions
1214-70	Esophageal Dilator Wall-mount Storage Cabinet	19"W x 4"D x 38"H
1214-71	Esophageal Dilator Horizontal Storage Cart	41"W x 23"D x 36"H
1214-72	Esophageal Dilator Vertical Storage Cart	24"W x 44"D x 24"H



Wall-Mount Storage Cabinet
#1214-70



Horizontal Storage Cart
#1214-71



Vertical Storage Cart
#1214-72

Replacement Parts & Accessories

Part No.	Description	Qty/Case
1214-73	HEPA Filter	1
1214-74	Storage Cabinet Cleaning Brush	10
1214-75	Soft Grip Gasket	1
1214-76	Fan	1
1214-77	Power Supply	1
1214-78	LED Light Strip, Single Side	1
1214-79	Swing Bar	1
1214-80	Dilator Storage Tray. Stores up to 12 bougies or dilators.	1


Hemostasis

Hemostasis

RePlay® Rotatable Hemostasis Clips

Single-use, Sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1170-01	11 mm width	230	2.8	5
1170-02	16 mm width	230	2.8	5
1170-03	18 mm width	230	2.8	5



RePlay® Rotatable Hemostasis Clip
#1170-01 thru 1170-03

Polypectomy & Tissue Acquisition

Biopsy Forceps

GIJaw™ Biopsy Forceps

Single-use, Sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1190-01	Serrated with Needle	230	2.8	10
1190-02	Serrated without Needle	230	2.8	10
1190-04	Jumbo, Serrated with Needle	230	3.2	10
1190-10	Smooth without Needle	160	2	10




Polypectomy Snares

Lasso™ Rotatable Polypectomy Snares

Single-use, Sterile

Part No.	Description	Indication	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1180-01	Oval 10 mm	Hot/Cold	230	2.8	10
1180-02	Oval 15 mm	Hot/Cold	230	2.8	10
1180-03	Oval 25 mm	Hot/Cold	230	2.8	10
1180-04	Oval 32 mm	Hot/Cold	230	2.8	10
1180-10	Hex 15 mm	Hot/Cold	230	2.8	10
1180-11	3-in-1 Hex 10/15/25 mm	Hot/Cold	230	2.8	10
1180-12	Oval 10 mm	Cold	230	2.8	10
1180-13	Oval 15 mm	Cold	230	2.8	10
1180-14	Hex 10 mm	Cold	230	2.8	10
1180-15	Hex 15 mm	Cold	230	2.8	10
4525	Universal Active Cord	—	—	—	1



Foreign Body Retrieval

Retrieval Nets

M-Catch® Retrieval Nets

Single-use, Non-sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1175-00	Standard Oval, 3 x 6 cm	160	2.8	5
1175-01	Standard Oval, 3 x 6 cm	230	2.8	5
1175-02	Standard Oval, 3 x 6 cm	350	2.8	5
1175-10	Polyp Oval, 3 x 6 cm	230	2.8	5




Grasping Forceps

EndoGrip™ Rotatable Grasping Forceps

Single-use, Non-sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1177-01	Rat Tooth-Alligator Combo	230	2.8	5



Endoscope Care & Accessories

Disposable Endoscope Valves & Accessories

PureFlo® Endoscope Valves & Accessories

Single-use

Part No.	Description	Piece/Kit	Compatibility	Sterility	Qty/Case
1195-21	Auxiliary Water Jet Connector*		Olympus	Sterile	100



PureFlo® Auxiliary Water Jet Connector
#1195-21

BullDog® Endoscope Valves & Accessories

Single-use

Part No.	Description	Piece/Kit	Compatibility	Sterility	Qty/Case
1195-22	Auxiliary Water Jet Connector*		Olympus	Sterile	100
1203-00	Biopsy Valve		Olympus/Fujinon	Non-sterile	100
1203-05	Biopsy Valve, bulk		Olympus/Fujinon	Non-sterile	200
1203-20	Biopsy Valve		Pentax	Non-sterile	100
1203-50	Irrigating Adapter			Non-sterile	50
1203-60	Air/Water Channel Cleaning Adapter		Olympus	Non-sterile	100
1203-70	Air/Water Valve, Suction Valve	2-Piece Kit	Olympus	Sterile	50
1203-71	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Olympus	Sterile	50
1203-72	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector*, Biopsy Valve	4-Piece Kit	Olympus	Sterile	50
1203-73	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector*, Biopsy Valve, Air/Water Cleaning Adapter	5-Piece Kit	Olympus	Sterile	40
1203-81	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Pentax	Sterile	50
1203-82	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector, Biopsy Valve	4-Piece Kit	Pentax	Sterile	50
1203-91	Air/Water Valve, Suction Valve, Biopsy Valve	3-Piece Kit	Fujinon	Sterile	50
1203-92	Air/Water Valve, Suction Valve, Auxiliary Water Jet Connector, Biopsy Valve	4-Piece Kit	Fujinon	Sterile	50



Olympus® 5-Piece Kit
#1203-73



Pentax® 4-Piece Kit
#1203-82

* Not compatible with Olympus® 140 Series GI Endoscopes.

For Olympus® 140/160/180/190/240/260/290 Series GI Endoscopes. Olympus is a registered trademark of Olympus Optical, Co., Ltd.

For Fujifilm® 700 Series Endoscopes. Fujinon is a registered trademark of Fujinon, Inc.

For Pentax® 90/110/110c Series Endoscopes. Pentax is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha.

Polyp Trap

Capture™ Polyp Trap

Single-use, Non-sterile

Part No.	Description	Qty/Case
1202-01	4 Chamber	25
1202-02	2-Tray	20



Capture® Polyp Trap:
4 Chamber
#1202-01



Capture® Polyp Trap:
2-Tray
#1202-02

Bite Blocks

CareGuard® Bite Blocks

Single-use, Non-sterile

Part No.	Size	Color	Style	Strap Option	Qty/Case
1196-48	48 FR	Green	With dental rim	Blue Strap	100
1196-60	60 FR	Pink	With dental rim and lateral ports	Blue Strap	100
1196-60XL	60 FR XL	Purple	With dental rim and lateral ports	Blue Strap	100
1197-60	60 FR	Pink	With dental rim and lateral ports	No Strap	100
1197-60XL	60 FR XL	Purple	With dental rim and lateral ports	No Strap	100
1198-48	48 FR	Green	With dental rim	Soft Strap	100
1198-60	60 FR	Pink	With dental rim and lateral ports	Soft Strap	100
1198-60XL	60 FR XL	Purple	With dental rim and lateral ports	Soft Strap	100



Careguard® Bite Block, 60 FR
#1196-60, 1197-60, 1198-60

CleanFreak® Endoscope Cleaning Products

Cleaning Brushes

CleanFreak® Endoscope Cleaning Brushes

Single-use, Non-sterile

Part No.	Description	Bristle Diameter (mm)	Sheath Diameter (mm)	Channel Diameter Range (mm)	Length (in/cm)	Qty/Case
1200-03	Extra Long Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-04	Dual Ended Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-05	Combination Brush (Valve & Channel Brush)	5-7	1.7	2.0-6.0	90.5/230	50
1200-06	Valve Brush	10				50
1200-07	Heavy Duty Standard Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-08	Heavy Duty Dual Ended Brush	5-7	1.7	2.0-6.0	90.5/230	50
1200-09	Heavy Duty Combination Brush (Valve & Channel Brush)	5-7	1.7	2.0-6.0	90.5/230	50



CleanFreak®
Heavy Duty
Combination Brush
#1200-09

Sponges & Pre-Cleaning Kit

CleanFreak® Endoscope Cleaning Sponges

Single-use, Non-sterile

Part No.	Description	Style	Color	Qty/Case
1201-01	Endoscope Sponge	Textured	Blue	200
1201-02	Endoscope Sponge	Cylindrical	Yellow	50
1201-03	Dry Imbibed Endoscope Sponge	Cylindrical	Yellow	50

CleanFreak® Pre-Cleaning Kit

Single-use, Non-sterile

Part No.	Description	Qty/Case
1201-20	Dry Imbibed Sponge with Water Basin	24



CleanFreak®
Pre-Cleaning Kit
#1201-20

Transport Pad

CleanFreak® Transport Pad

Single-use, Non-sterile

Part No.	Description	Size	Qty/Case
1201-90	Transport Pad	40" x 29"	25



CleanFreak®
Transport Pad
#1201-90

Endoscopy Kits

Basic Endoscopy Kits

Diversatek Healthcare will customize a kit to your exact specifications. Please contact your Sales Representative for details on customizing your own kit.



Basic Endoscopy Kit

Lighted Bougies

InnerVision® Transillumination System (Lighted Bougie)

InnerVision® Fiber Optic Cable

Reusable, Non-sterile

Part No.	Size	Compatible Connector	Qty/Case
1450-40	40 Fr	Circon-ACMI/Olympus	1
1450-50	50 Fr	Circon-ACMI/Olympus	1
1450-56	56 Fr	Circon-ACMI/Olympus	1
1455-40	40 Fr	Storz	1
1455-50	50 Fr	Storz	1
1455-56	56 Fr	Storz	1

InnerVision® Fiber Optic Tips

Single-use, Sterile

Part No.	Size	Qty/Case
1465-40	40 Fr, Blunt	5
1465-50	50 Fr, Tapered	5
1465-56	56 Fr, Tapered	5

Detachable InnerVision® Tips must be used with the corresponding InnerVision® Cable

Wolf/Dyonics Adapter for Circon-ACMI/Olympus Light Source

Reusable, Non-sterile

Part No.	Qty/Case
1450-02	1

Wolf/Dyonics Adapter, Storz Connector and Circon-ACMI/Olympus Connector.