


# CareGuard® Bite Blocks

## Protection and Confidence

Color-coded bite blocks allow for efficient size identification,  
and the single-use design supports patient safety.



## Bite Blocks

### Reliable Performance

A sturdy dental retention rim helps the bite block maintain its position to protect both the patient and the endoscope.



#### 48 FR Size

Identified by a distinct green color, this bite block provides enhanced comfort for smaller-sized patients.

#### 60 FR Size

A vibrant pink color provides for easy identification during procedure setup with this standard-sized bite block.

#### 60 FR XL Size

The largest bite block features enlarged lateral ports to simplify scope positioning and improve access for additional suctioning equipment.

#### Adjustable Options

Single-use bite blocks are available with easy-to-use adjustable straps.



**Diversatek™**  
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: [orders@diversatekhc.com](mailto:orders@diversatekhc.com)

General Inquiries: [sales@diversatekhc.com](mailto:sales@diversatekhc.com)

Global Corporate Office:  
102 East Keefe Avenue  
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://diversatekhealthcare.com)   


©2025 Diversatek, Inc. | All Rights Reserved. 080-0032 12/25

Access On-Demand  
Education



[diversateku.com](https://diversateku.com)

Scan for  
Product Details



[divhc.co/bite](https://divhc.co/bite)