


Esophageal Dilator Storage

Superior Storage

Innovative storage solutions keep esophageal dilators fully enclosed and free from contact with one another.



Dilator Storage Solutions

Protect and Preserve

Safeguard the integrity of esophageal dilators and support infection prevention with purposefully designed storage systems.

405 nm LED Lights

Help to deactivate bacteria and maintain tip integrity with LED lighting.


HEPA Filters & Fan

Minimize exposure to particulates and airborne contaminants with forced air circulation and a positive pressure environment.



Complete Portfolio

Three distinct storage solutions accommodate your needs to maintain clean, straight, and secure esophageal dilators.



Wall-Mount Cabinet



Horizontal Cart



Vertical Cart

Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://diversatekhealthcare.com)   

©2025 Diversatek, Inc. | All Rights Reserved. 080-0077 12/25

Access On-Demand
Education



diversateku.com

Scan for
Product Details



divhc.co/storage