


Lasso™ Rotatable Polypectomy Snares

Optimal Control


Achieve precise placement with responsive rotation and a balanced wire design to confidently secure and remove polyps.



Polypectomy Snares

Secure Polyp Capture

The braided wire design securely holds tissue and delivers a clean, precise cut across cold and hot/cold snares.



Thin Braided Wire

Rotatable **cold-only** snare with a thin (0.024") wire ensures a clean cut for flat and sessile polyps, enabling more complete resections.

Full Portfolio

A comprehensive selection of sizes and shapes for both **hot and cold** applications.



Full 360° Rotation

An ergonomic handle and short throw enable controlled actuation with 1:1 rotation.



3-in-1 Snare

Utilize a single snare to accommodate multiple polyp sizes. This hot/cold hex snare is adjustable to three distinct diameters, reducing the need for multiple snares in a procedure.

Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://diversatekhealthcare.com)   

©2025 Diversatek, Inc. | All Rights Reserved. 080-0076 12/25

Access On-Demand
Education



diversateku.com

Scan for
Product Details



divhc.co/snares