


M-Catch® Retrieval Nets

Reliable Retrieval

Interwoven netting reinforced around the wire loop delivers strength and security to capture foreign bodies and polyps.




Retrieval Nets

Precise Control

Designed for smooth actuation, the handle and netting enable accurate navigation, positioning, and secure retrieval.

Netting Options

Utilize the polyp net with larger openings or the standard net with a tighter weave to meet retrieval needs.



Interwoven Netting

Robust netting allows repeated opening and closing without compromising a secure, full closure around various foreign bodies and polyps.



PTFE Tubing

Durable construction resists kinking and bending to enable smooth insertion into the endoscope channel.



Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://diversatekhealthcare.com)   


©2025 Diversatek, Inc. | All Rights Reserved. 080-0072 12/25

Access On-Demand
Education



diversateku.com

Scan for
Product Details



divhc.co/nets