

Viper® 3-Stage Balloon Dilator

Superior Strength

A balloon engineered to deliver enhanced radial force and uniform expansion across all three stages.



Balloon Dilators

Visible Precision

Visually clear balloon material with high-contrast marker bands enables precise placement and controlled dilation.

PET Balloon Technology

Provides optimal radial strength and uniform dilation with a sustained diameter — no repeat inflations needed.

Nitinol Core Wire

Kink-resistant nitinol core enhances pushability and maintains shape for precise control with the fixed wire balloon dilator.

Comprehensive Offering

Fixed and wire-guided options available in a full range of sizes (6 mm to 20 mm) from a trusted source.

Atraumatic Tip

Soft, flexible tips are designed for atraumatic access to minimize the risk of mucosal injury.

Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://www.DiversatekHealthcare.com)   

©2025 Diversatek, Inc. | All Rights Reserved. 080-0086 12/25

Access On-Demand
Education



[diversateku.com](https://www.diversateku.com)

Scan for
Product Details



divhc.co/viper