

InnerVision® Transillumination System

Precise Identification & Visualization

Flexible, lighted bougie with a clear silicone tip emits low-heat, radial light for consistent illumination during surgical procedures.



Lighted Bougie

Safe Navigation

Silicone construction provides flexibility, safety, and radial illumination in a detachable, disposable tip.



Tip Configuration

Available in blunt and tapered tips to accommodate varying surgical needs and preferences.

Multiple Sizes

Single-use sterile tips come in 40 FR, 50 FR, and 56 FR measurements.

Compatible Connectors

Engineered with a robust snap-fit design that securely connects the tip to fiber optic cables using a 250–300 watt light source.



Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://www.DiversatekHealthcare.com)   

©2025 Diversatek, Inc. | All Rights Reserved. 080-0011 12/25

Access On-Demand
Education



[diversateku.com](https://www.diversateku.com)

Scan for
Product Details



divhc.co/light