

EZ-Ject™ Injection Needle

Dependable Delivery

Designed with a double-beveled, tapered needle and one-step locking mechanism for reliable delivery of viscous solutions.



Endoscopic Injection Needle

Versatile Needle for Everyday Injections

The injection needle is designed to deliver a range of viscous solutions for lift-assisted polypectomy, hemostasis, and endoscopic tattooing.



Clear Catheter Design

A clear catheter with an optimized inner diameter enables uninterrupted flow of viscous solutions and allows for real-time visualization during injection.



One-Step Lock

The contoured handle features a strategically placed thumb-and-forefinger groove that enables safe, easy needle actuation. An audible and tactile clicking mechanism provides confident confirmation that the needle is securely in place.



Double-Beveled Needle

The 5 mm needle features a tapered tip for smooth submucosal puncture with minimal force, reducing tissue trauma. Available in 22 and 25 gauge.

Diversatek™
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com

General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

[DiversatekHealthcare.com](https://diversatekhealthcare.com)   

©2026 Diversatek, Inc. | All Rights Reserved. 080-0273 01/26

Access On-Demand
Education



diversateku.com

Scan for
Product Details



divhc.co/inject