

ZepHr® Reflux Monitoring System

Complete Assessment

Capture acid and non-acid reflux episodes along with impedance, providing a deeper understanding of patient symptoms to guide more informed treatment decisions.



Reflux Monitoring

Actionable Insights


The ZepHr® Reflux Monitoring System captures 24-hour pH and impedance data from the lower esophagus to the laryngopharynx. Aligned with the latest Lyon Consensus normal values, ZepHr delivers clinically meaningful insights to support accurate evaluation and confident treatment decisions.

ComforTEC® Z/pH & pH Probes

Featuring a comprehensive selection of 29 probe configurations, including:

- Z/pH including laryngopharyngeal (LPR) (impedance and acid)
- pH only (acid)

Each probe is constructed from specially formulated polyurethane that softens at body temperature to maximize patient comfort.



ZepHr® Recorder

- **High-contrast**, backlit display makes it easy-to-read
- **Intuitive buttons** simplify marking meals, sleep, and diary entries
- **Compact** and light for enhanced patient comfort and compliance

Zvu® Functional GI Software

Zvu Software with AutoSCAN® technology automatically analyzes pH and impedance data to identify reflux events, classify acid and non-acid events on or off PPIs, track bolus movements, and calculate key metrics such as **acid exposure time**, **MNBI**, and the **DeMeester score**. Customizable reports and clear visualization tools further enhance efficiency and diagnostic confidence.



Diversatek
Healthcare

+1 800 558 6408 or +1 414 265 7620

Orders: orders@diversatekhc.com
General Inquiries: sales@diversatekhc.com

Global Corporate Office:
102 East Keefe Avenue
Milwaukee, WI 53212 USA

DiversatekHealthcare.com   

©2026 Diversatek, Inc. | All Rights Reserved. 080-0080 02/26

Access On-Demand
Education



diversateku.com

Scan for
Product Details



divhc.co/zephr